





BinaxNOW®
RSV Card

Intended Use

The Alere BinaxNOW® RSV Card is a rapid immunochromatographic assay for the qualitative detection of respiratory syncytial virus (RSV) fusion protein antigen in nasal wash and nasopharyngeal swab specimens from symptomatic patients. This test is intended for *in vitro* diagnostic use to aid in the diagnosis of respiratory syncytial virus infections in neonatal and pediatric patients under the age of 5. It is recommended that negative test results be confirmed by cell culture.

Summary and Explanation of the Test

RSV is a common cause of upper and lower respiratory tract infections and the major cause of bronchiolitis and pneumonia in infants and children. Infections and outbreaks due to RSV typically occur yearly in the fall, winter and spring. While RSV can cause significant respiratory illness in older children and adults, the disease tends to be milder in these populations than in infants and young children.

Rapid identification and diagnosis of RSV has become more important due to the availability of effective antimicrobial therapy. Rapid identification can lead to reduced hospital stays, reduction in antimicrobial use, and reduction in the cost of hospital care.¹

The Alere BinaxNOW® RSV Card provides a simple, rapid method for the diagnosis of RSV using nasal wash and nasopharyngeal swab specimens. The easy-to-use format and rapid results allow for its use in “STAT” testing where it can provide invaluable information to assist with treatment and hospitalization decisions.

Principles of the Procedure

The Alere BinaxNOW® RSV Card is an immunochromatographic membrane assay used to detect RSV fusion protein antigen in nasal wash and nasopharyngeal swab specimens. Anti-RSV antibody, the Sample Line, is adsorbed onto nitrocellulose membrane. Control antibody is adsorbed onto the same membrane as a second stripe. Both anti-RSV and control antibodies are conjugated to visualizing particles that are dried onto an inert fibrous support. The resulting conjugate pad and the striped membrane are combined to construct the test strip. This test strip is mounted on the right side of a cardboard, book-shaped hinged test card.

Swab samples (controls and patients) require a preparation step, in which the sample is eluted off the swab into an appropriate solution. Nasal wash samples do not require any preparation.

To perform the test, the sample to be tested is added to the white pad at the top of the test strip, and the test card is closed. RSV antigen present in the sample reacts to bind anti-RSV conjugated antibody. The resulting antigen-conjugate complexes are captured by immobilized anti-RSV antibody, forming the Sample Line. Immobilized Control Line antibody captures a visualizing conjugate, forming a pink Control Line. The Control Line is blue in a card that has not been tested.

Test results are interpreted by the presence or absence of visually detectable pink-to-purple colored lines. A positive test result, read at 15 minutes, will include the detection of both a Sample Line and a Control Line. A negative test result, read at 15 minutes, will produce only a Control Line, indicating that RSV antigen was not detected in the sample. Failure of the Control Line to appear, or the Control Line remaining blue, indicates an invalid assay, whether the Sample Line is present or not.

Reagents and Materials

Materials Provided

Refer to illustrations on pull-out flap.

- 1 **Test Cards:** A membrane coated with mouse antibody specific for RSV antigen and with control line antibody is combined with mouse anti-RSV and control line antibody conjugates in a hinged test card. The membrane of an untested card contains a blue line at the control line area.
- 2 **Transfer Pipettes:** Fixed volume (100 µl) transfer pipettes used to transfer sample to the test cards. Use only pipettes provided by Alere or a calibrated pipette capable of delivering 100 µl sample volume.
- 3 **Positive Control Swab:** Inactivated RSV dried onto swab.
- 4 **Negative Control Swab:** Inactivated *Streptococcus* Group A dried onto swab.
- 5 **Elution Solution Vials for Control Swabs:** Vials containing elution solution used to prepare the Control Swabs. Do not use other elution solutions with the Alere BinaxNOW® RSV Card. 

Materials Not Provided

Clock, timer or stopwatch; nasal wash collection containers, nasopharyngeal swabs, transport media

Optional Item

Nasopharyngeal (Np) Swab Specimen Accessory Pack:

- 6 **Nasopharyngeal Swabs:** Sterile swabs for use in the Alere BinaxNOW® RSV Card.
- 7 **Elution Solution Vials for Swab Specimens:** Vials containing elution solution used to prepare swab specimens for testing. 

Precautions

1. For *in vitro* diagnostic use.
2. Leave test card sealed in its foil pouch until just before use.
3. Do not use kit past its expiration date.
4. Do not mix components from different kit lots.
5. The white sample pad at the top of the test strip contains reagents that extract the target antigen from the virus. To ensure optimum performance, add the sample **SLOWLY** to the **MIDDLE** of this pad such that all of the sample volume absorbs into the pad.
6. The RSV Positive Control Swab has been prepared from RSV-infected tissue culture cells that have been inactivated and subsequently tested by bioassay procedures. Use universal precautions when performing the assay. Samples may be infectious. Proper handling and disposal methods should be established according to local, state, and federal regulations.
7. **INVALID RESULTS** can occur when an insufficient volume of specimen is added to the test card. To ensure delivery of an adequate volume, make certain that the lower shaft of the transfer pipette is full and does not contain air spaces before dispensing contents of the pipette onto the Sample Pad of the card. If air spaces are present, expel the specimen back into the container by squeezing the top bulb and redraw the specimen into the pipette. Use a new pipette if necessary.
8. When testing nasal wash samples, avoid viscous areas of the sample when drawing specimen into the transfer pipette. If the pipette becomes clogged, such that the lower shaft of the pipette is not full, expel the specimen back into container by squeezing the top bulb and redraw the specimen into the pipette. Use a new pipette if necessary.

9. Polyester, rayon, foam, cotton or flocked flexible shaft nasopharyngeal swabs have been evaluated and found to be acceptable for use in the Alere BinaxNOW® RSV Card. Do not use calcium alginate nasopharyngeal swabs in the Alere BinaxNOW® RSV Card.
10. All transfer pipettes and elution vials are single use items – do not use with multiple specimens.
11. Elution solution contains Triton® X-100 and ProClin® 300. Warning; may cause an allergic skin reaction, causes serious eye irritation. ⚠
12. Safety Data Sheets for this product are available upon request.
13. Follow your national, regional, and local ordinances accordingly for waste disposal regulations.

Storage and Stability

Store kit at 2-30°C. The Alere BinaxNOW® RSV Card kit and reagents are stable until the expiration date printed on kit box.

Specimen Collection and Handling

Use fresh nasopharyngeal swab or nasal wash specimens for optimal Alere BinaxNOW® RSV Card performance.

Nasopharyngeal Swabs:

Polyester, rayon, foam, cotton and flocked nasopharyngeal swabs, all with flexible shafts, have been evaluated and found to be acceptable for use in the Alere BinaxNOW® RSV Card. Add swab specimens to 0.5 - 3.0 ml of a suitable liquid transport system within one hour of collection. If immediate testing is not possible, eluted liquid swab samples can be stored at room temperature for up to 4 hours, or at 2-8°C for up to 48 hours, before testing. Allow samples to warm to room temperature and swirl gently before testing.

Nasal Washes:

Collect wash samples in standard collection cups. Use procedures appropriate for the age of the patient. If immediate testing is not possible, washes can be stored at room temperature for up to 4 hours, or at 2 - 8°C for up to 24 hours, before testing. Allow samples to warm to room temperature and swirl gently before testing.

Wash samples can be placed in up to 3 ml of a suitable liquid transport system prior to testing in the Alere BinaxNOW® RSV Card. Use of transport media will result in dilution of wash samples. This dilution may lower the overall test sensitivity as reported in the Performance Section.

Transport Media:

The following transport media were tested and found to be acceptable for use in the Alere BinaxNOW® RSV Card.

Amies Media	Phosphate Buffer Solution
Brain Heart Infusion Broth	Saline
Dulbecco Medium	Stuart's Media
Hank's Balanced Salt Solution	Tryptose Phosphate Broth
M4 Media	UTM-RT Media
M4-RT Media	Veal Infusion Broth
M5 Media	

Quality Control

Daily Quality Control:

The Alere BinaxNOW® RSV Card has built-in procedural controls. For daily quality control, Binax suggests that you record these controls for each test run.

Procedural Controls

- An untested card has a blue line at the "Control" position. If the test flows and the reagents work, this blue line will always turn pink in a tested card.
- The clearing of background color from the result window is a negative background control. The background color in the window should be light pink to white within 15 minutes. Background color should not hinder reading of the test.

External Positive and Negative Controls:

Good laboratory practice suggests the use of positive and negative controls to ensure that:

- test reagents are working, and
- the test is correctly performed.

Alere BinaxNOW® RSV Card kits contain Positive and Negative Control Swabs. These swabs will monitor for substantial reagent failure. The Positive Control will not ensure precision at the assay cut-off. Test these swabs once with each new shipment received. Other controls may be tested in order to conform with:

- local, state and/or federal regulations,
- accrediting groups, and/or,
- your lab's standard Quality Control procedures.

Refer to 42 CFR 493.1256 for guidance on proper QC practices (U.S. customers only).

If the correct control results are not obtained, do not report patient results. Contact Technical Service during normal business hours.

Sample Preparation Procedure

Nasal Washes:


Nasal wash samples do not need preparation. Go to Test Procedure.

Nasopharyngeal Swabs:

Elute swab in 0.5-3.0 ml of a liquid transport system by rotating swab in the liquid. Go to Test Procedure. Note: If eluting swab in the Alere BinaxNOW® Elution Solution, follow the procedure for Control Swab preparation below.


Control Swabs:

1. The test kit contains test vials pre-filled with elution solution. Twist off the test vial cap.
2. Put swab to be tested into test vial. Rotate the swab three (3) times in the liquid.
3. Press the swab against the side of the vial and turn as you remove it from the vial. This removes sample from the swab.
4. Discard the swab.
5. Test the liquid sample (from the test vial) in the Alere BinaxNOW® RSV Card as soon as possible. Go to Test Procedure.

Test Procedure

1. Remove card from the pouch just prior to testing and lay flat on work bench.
2. Fill pipette by firmly squeezing the top bulb and placing pipette tip into liquid sample. Release bulb while tip is still in sample. This will pull liquid into the pipette. Make sure there are no air spaces in the lower part of the pipette.
3. See arrow on test card to find **White** Sample Pad. **SLOWLY** add entire contents (100 μ l) of pipette to the **MIDDLE** of this pad by squeezing the top bulb. **DO NOT** add sample to the pink/purple colored pad.
4. Immediately peel off adhesive liner from the test card. Close and securely seal the card. Read result in window 15 minutes after closing the card. Results read before or after 15 minutes may be inaccurate.

Note: When reading test results, tilt the card to reduce glare on the result window if necessary.

Result Interpretation


For a **NEGATIVE SAMPLE**, the BLUE Control Line in the lower half of the window turns a PINK-TO-PURPLE color. No other line appears.

For a **POSITIVE SAMPLE**, the BLUE Control Line turns a PINK-TO-PURPLE color. A second PINK-TO-PURPLE Sample Line appears above it.

A test is **INVALID** if the Control Line remains blue or is not present at all. Repeat invalid tests with a new test card. Call Technical Service if the problem persists.

Reporting of Results

Result Suggested Report

- Positive** Positive for RSV antigen. A positive result may occur in the absence of viable virus.
- Negative** Negative for RSV antigen. Infection due to RSV cannot be ruled out. The antigen in the sample may be below the detection limit of the test. Alere suggests culture of negative samples.

Limitations

The Alere BinaxNOW[®] RSV Card detects both viable and non-viable RSV. Test performance depends on antigen load in the specimen and may not correlate with cell culture performed on the same specimen.

Inadequate specimen collection or low levels of virus shedding may result in suboptimal performance and may yield false negative results.

A negative test result does not exclude infection with RSV nor is it intended to rule out other microbial-caused respiratory infections. Therefore, the results obtained with the Alere BinaxNOW[®] RSV Card should be used in conjunction with clinical findings to make an accurate diagnosis.

Alere BinaxNOW[®] RSV Card performance has not been evaluated in patients who have been treated with palivizumab. However, an analytical study has demonstrated that palivizumab interferes with the ability of the Alere BinaxNOW[®] RSV Card to detect RSV.

The potential for interference from anti-microbials and interferon has not been established.

Monoclonal antibodies may not detect all antigenic variants or new strains of RSV.

Expected Values

The prevalence of RSV varies from year to year, with outbreaks typically occurring during the fall and winter months. The rate of positivity found in RSV testing is dependent on many factors including the method of specimen collection, the test method used, geographic location, and the disease prevalence in specific localities. In the 2002 Alere[™] clinical study, the average prevalence of RSV was 2% in wash samples and 4% in nasopharyngeal swab samples. Prevalence of RSV in nasopharyngeal swab samples collected during the 2003 Alere[™] clinical study was 21%.

Performance Data

ANALYTIC STUDIES Analytical Reactivity:

There are 2 known subgroups of respiratory syncytial virus (RSV) and both contain the conserved fusion protein targeted by the Alere BinaxNOW[®] RSV Card.² Six (6) subgroup A clinical isolates and five (5) subgroup B clinical isolates tested positive in the Alere BinaxNOW[®] RSV Card at concentrations ranging from 1.56×10^3 TCID₅₀/ml to 5.00×10^4 TCID₅₀/ml.

Note: The reported TCID₅₀/ml levels are dependent on a number of factors including the cell culture lines used, the number of passages performed and the efficiencies of the various isolates.

Analytical Specificity (Cross-Reactivity):

To determine the analytical specificity of the Alere BinaxNOW[®] RSV Card, 48 commensal and pathogenic microorganisms (28 bacteria and 20 viruses) that may be present in the nasal cavity or nasopharynx were tested. All of the following microorganisms were negative when tested at concentrations greater than 1×10^5 TCID₅₀/ml (viruses) or greater than 1×10^8 organisms/ml (bacteria). Metapneumovirus was tested at 2×10^3 TCID₅₀/ml and did not cross-react.

Bacteria	Virus
<i>Acinetobacter</i>	Adenovirus 5*
<i>Bordetella pertussis</i>	Adenovirus 7*
<i>Candida albicans</i>	CMV*
<i>Enterococcus faecalis</i>	Coronavirus*
<i>Escherichia coli</i>	Coxsackie B4*
<i>Gardnerella vaginalis</i>	Influenza A 2/Japan/305/57
<i>Haemophilus influenzae</i>	Influenza A/Hong Kong/8/68
<i>Klebsiella pneumoniae</i>	Influenza A/Aichi/68
<i>Lactobacillus casei</i>	Influenza A/PR/8/34
<i>Legionella pneumophila</i>	Influenza A/Victoria/3/75
<i>Listeria monocytogenes</i>	Influenza A1/FM/1/47
<i>Moraxella catarrhalis</i>	Influenza B Allen/45
<i>Neisseria gonorrhoeae</i>	Influenza B Lee/40
<i>Neisseria meningitidis</i>	Influenza B Mass/3/66
<i>Neisseria sicca</i>	Influenza B Maryland/1/59
<i>Neisseria subflava</i>	Influenza B Taiwan/2/62
<i>Proteus vulgaris</i>	Metapneumovirus
<i>Pseudomonas aeruginosa</i>	Parainfluenza 1*

Bacteria	Virus
<i>Serratia marcescens</i>	Parainfluenza 2*
<i>Staphylococcus aureus</i>	Parainfluenza 3*
<i>Staphylococcus aureus</i> (Cowan protein A producing strain)	
<i>Staphylococcus epidermidis</i>	
<i>Streptococcus</i> , Group A	
<i>Streptococcus</i> , Group B	
<i>Streptococcus</i> , Group C	
<i>Streptococcus</i> , Group F	
<i>Streptococcus mutans</i>	
<i>Streptococcus pneumoniae</i>	

*These viral strains were obtained from American Type Culture Collection (ATCC) with titer information, and the titers were not verified by Alere.

Interfering Substances:

The following substances, naturally present in respiratory specimens or that may be artificially introduced into the nasal cavity or nasopharynx, were evaluated in the Alere BinaxNOW® RSV Card at the concentrations listed and were found not to affect test performance.

Substance	Concentration
Whole blood	2%
3 OTC mouthwashes	25%
3 OTC throat drops	25%
3 OTC nasal sprays	25%
4-acetamidophenol	10 mg/ml
Acetylsalicylic acid	20 mg/ml
Albuterol	20 mg/ml
Chlorpheniramine	5 mg/ml
Dextromethorphan	10 mg/ml
Diphenhydramine	5 mg/ml
Guaiaicol glycerol ether	20 mg/ml
Oxymetazoline	10 mg/ml
Phenylephrine	100 mg/ml
Phenylpropanolamine	20 mg/ml
Rebetol®	500 ng/ml
Relenza®	20 mg/ml
Rimantadine	500 ng/ml
Tamiflu®	100 mg/ml

CLINICAL STUDIES

Nasal Wash – Clinical Specificity (Prospective Study):

The performance of the Alere BinaxNOW® RSV Card was compared to cell culture in a multi-center study conducted during the 2002 Flu season at physician offices and clinics located throughout the United States. Nasal wash specimens were collected from children and adults presenting with RSV-like symptoms for 3 days or less and evaluated in the Alere BinaxNOW® RSV Card. The population tested was approximately 46% female and 54% male. Patients were not included in the study if they had received an influenza vaccine within 6 months of the enrollment period, or if they had taken either an influenza or RSV medication within 30 days of the enrollment period. There were no invalid tests reported.

One hundred ninety-one (191) nasal wash specimens were tested at 4 different test sites. Alere BinaxNOW® RSV Card overall specificity was 98%, and overall test agreement was 98%. Ninety-five percent (95%) confidence intervals are listed below.

		Wash Viral Culture	
		+	-
Alere BinaxNOW® RSV Card	+	3	3
Result	-	1	184
		95% CI	
Specificity	=	98% (184/187)	(95.4% - 99.4%)
Overall Agreement	=	98% (187/191)	(94.7% - 99.1%)

The Alere BinaxNOW® RSV Card performed similarly at the 4 test sites as shown in the table below.

	Positive Points	Specificity		
		Alere BinaxNOW® RSV Card / Culture	%	95% CI
Site I	1/1	91/94	97	91.1 - 98.8
Site II	2/3	83/83	100	95.7 - 100
Site III	0/0	6/6	100	59.0 - 99.6
Site IV	0/0	4/4	100	47.8 - 99.5

Nasal Wash – Clinical Sensitivity and Specificity (Retrospective Study):

Since there were a low number of positive culture confirmed RSV infections generated during the prospective study, a retrospective study was conducted as follows. Nasal wash specimens from 47 viral culture positive RSV patients and 12 viral culture negative RSV patients were evaluated in the Alere BinaxNOW® RSV Card. All of the samples were obtained from a large university medical center and had been collected from patients living in the northeastern region of the US. The population tested was approximately 49% male and 51% female.

Alere BinaxNOW® RSV Card sensitivity was 89%, while test specificity was 100%. Overall test agreement was 92%. Ninety-five percent (95%) confidence intervals are listed below.

		Wash Viral Culture	
		+	-
Alere BinaxNOW® RSV Card	+	42	0
Result	-	5	12
		95% CI	
Sensitivity	=	89% (42/47)	(77.3% - 95.3%)
Specificity	=	100% (12/12)	(75.3% - 99.8%)
Overall Agreement	=	92% (54/59)	(81.6% - 96.2%)

Nasopharyngeal Swab - Sensitivity and Specificity (Prospective Study):

The performance of the Alere BinaxNOW® RSV Card on nasopharyngeal swab specimens was compared to cell culture/DFA in a multi-center US study conducted during the 2002 and 2003 Flu seasons. Nasopharyngeal swab specimens were collected from children presenting with RSV or flu-like symptoms. All swab samples were placed in 0.5-3 ml of viral transport media prior to evaluation in the Alere BinaxNOW® RSV Card. The population tested was 43% female and 57% male.

One hundred and seventy-nine (179) nasopharyngeal swab specimens were tested. There were no invalid tests reported. Alere BinaxNOW® RSV Card sensitivity, specificity and overall agreement as compared to culture/DFA were 93%. Ninety-five percent (95%) confidence intervals are listed below.

Nasopharyngeal Swab

Culture / DFA

Alere BinaxNOW® RSV Card Result		Culture / DFA	
		+	-
+		25	10
-		2	142

95% CI

Sensitivity	=	93% (25/27)	(76.5% - 97.7%)
Specificity	=	93% (142/152)	(88.3% - 96.4%)
Overall Agreement	=	93% (167/179)	(88.6% - 96.1%)

The Alere BinaxNOW® RSV Card performed similarly at all test sites as shown in the tables below.

SENSITIVITY				SPECIFICITY			
Site	#	%	95% CI	Site	#	%	95% CI
1	14/15	93	69.8-98.4	1	69/74	93	85.1-97.0
2	9/10	90	58.7-97.7	2	20/23	87	67.6-95.3
3	0/0	NA	NA	3	16/18	89	66.9-96.6
4	2/2	100	29.2-99.2	4	37/37	100	90.7-99.9

Reproducibility Study:

A blind study of the Alere BinaxNOW® RSV Card was conducted at 3 separate sites using panels of blind coded specimens containing negative, low positive, and moderate positive samples. Participants tested each sample multiple times on 3 different days. One hundred percent (100%) of the 234 samples tested produced the expected result yielding a 95% confidence interval of 98.4 - 100%.

Ordering and Contact Information

Reorder numbers:

#430-000:	Alere BinaxNOW® RSV Card (42 test kit)
#430-022:	Alere BinaxNOW® RSV Card (22 test kit)
#400-065:	Alere BinaxNOW® Nasopharyngeal Swab Accessory Pack (20 swab kit)
#430-080	Alere BinaxNOW® RSV Control Swab Kit

US 1 877 441 7440
OUS +1-321-441-7200

Technical Support Advice Line

Further information can be obtained from your distributor, or by contacting Alere™ Technical Support on:

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© 2015 Alere. All rights reserved.
The Alere Logo, Alere and BinaxNOW are trademarks of the Alere group of companies.
All trademarks referenced are trademarks of their respective owners.
**Alere Scarborough, Inc.**