



Liquid with a surprise

HomeLox makes and stores liquid oxygen at home

PHILIPS

RESPIRONICS

sense **and** simplicity

Liquid in a whole



new light



The innovative HomeLox portable liquid oxygen system from Philips Respironics enables your patients to now generate and store liquid oxygen right in their own home. This groundbreaking system offers the long-lasting portability of liquid oxygen in a way that greatly reduces the cost and complexity of traditional systems.

Think liquid

Whether you are currently committed to providing liquid oxygen to your patients, outsourcing liquid to other providers, or turning away liquid referrals altogether, the HomeLox system can benefit your bottom line.

Outsourcing liquid?

Gain control by providing liquid yourself—without the investment required in traditional systems.

Current liquid oxygen provider?

Transition away from the expensive equipment and recurrent delivery costs while providing patients the ability to generate oxygen in their homes.

Not offering liquid?

Grow your oxygen business by offering a unique liquid oxygen solution to clinicians. Maintain your relationships and stop turning away referrals by including liquid oxygen in your offerings.

Better for your bottom line

HomeLox eliminates many of the significant costs that are traditionally associated with delivered liquid oxygen including:

- Recurrent deliveries
- Specialized vehicles
- Testing equipment
- Storage containers
- Filling-related service calls

The easy-to-use system includes the HomeLox home liquid oxygen generation and storage unit and the highly portable, long-lasting GoLox portable device.



Let us fill you in

The HomeLox portable liquid oxygen system uses a new hands-free clean-fill process that reduces the freezing, overfilling, or under filling typically found in traditional liquid systems.

Our clean-fill process:

- Automatically stops filling once GoLox is full
- Maximizes portable duration
- Eliminates calls from patients who struggle with traditional liquid filling



GoLox portable liquid oxygen device

- Small and lightweight (3.8 lbs filled)
- No batteries, no maintenance required
- Lasts up to 10 hours (on a setting of 2 at 20 breaths per minute)
- Uses a single-lumen cannula



on innovation



Three steps to making liquid oxygen at home

Step 1: Concentrate

Using a high-quality compressor and the same sieve bed material found in the EverFlo stationary oxygen concentrator, HomeLox collects room air and converts it into prescription-grade oxygen.

Step 2: Refrigerate/store

The oxygen is then cooled to -279°F (-173°C) using a proprietary refrigerant and standard refrigeration technology typically found in industrial and household applications. The chilled oxygen is converted into a liquid and then stored.

Step 3: Refill

When it's time to refill GoLox, simply place the unit on top of HomeLox, lock the unit, and pull the filling lever. This hands-free process addresses the difficulties found with traditional liquid systems that require the use of one or both hands to fill. When full, filling stops automatically.

Easy to maintain

For providers, the HomeLox system is low-maintenance with only one filter change required every two years. For patients, HomeLox is maintenance free.

HomeLox specifications	
Size	35 3/8" h x 18 1/2" w x 18 1/2" d
Weight	109 lbs (empty), 116 lbs (full)
Liquid oxygen generation	3 liters per day (typical)
Liquid oxygen storage capacity	3 liters
GoLox fill time	< 1 minute
Normal evaporation time	.5 liters per day
Warranty	2 years or 6,000 hours (whichever comes first)
Reimbursement code	E0433 Portable liquid oxygen system, rental Home liquefier used to fill portable liquid oxygen containers, includes portable containers, regulator, flow meter, humidifier, cannula or mask and tubing, with or without supply reservoir and contents gauge

GoLox specifications		
Liquid oxygen capacity	.375 liters (0.9 lbs)	
Gaseous equivalent	308 liters	
Weight	3.8 lbs (filled), 2.9 lbs (empty)	
Dimensions	7.25" h x 7" w x 3.75" d	
Duration	10 hrs (typical setting 2, 20 BPM) 2.5 hrs continuous flow	
Normal evaporation rate	1.0 lb per day (typical), 1.5 lbs (max)	
Warranty	Device - 1 year Vacuum - 5 years	
Flow settings and amount	Delivery dosage settings	Volume (ml/breath)
	1	12
	1.5	18
	2	24
	2.5	30
	3	36
	4	48

Ordering information

Description	Part number
HomeLox portable liquid oxygen system (Includes the HomeLox home liquefaction unit and GoLox portable liquid oxygen device)	1057020
HomeLox home liquefaction unit	1047777
GoLox portable liquid oxygen device (Includes the HomeLox compatible GoLox, 4' single-lumen cannula, carry bag shoulder strap, and user manuals)	1019692

Accessories	Part number
Carry bag (A carry bag for the GoLox with a detachable shoulder strap that can be used as a shoulder or waist strap)	1018019
HomeLox user's manual	1047833
GoLox user's manual (HomeLox compatible version)	1052483

**Philips Healthcare is part of
Royal Philips Electronics**

How to reach us

www.philips.com/healthcare
healthcare@philips.com

Asia
+49 7031 463 2254

Europe, Middle East, Africa
+49 7031 463 2254

Latin America
+55 11 2125 0744

North America
+1 425 487 7000
800 285 5585 (toll free, US only)

Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668

Customer Service
+1 724 387 4000
800 345 6443 (toll free, US only)

Philips Respironics International
Headquarters
+33 1 47 28 30 82

Philips Respironics Asia Pacific
+65 6882 5282

Philips Respironics Australia
+61 (2) 9666 4444

Philips Respironics China
+86 021 24127311

Philips Respironics Deutschland
+49 8152 93 06 0

Philips Respironics France
+33 2 51 89 36 00

Philips Respironics Italy
+39 039 203 1

Philips Respironics Sweden
+46 8 120 45 900

Philips Respironics Switzerland
+41 6 27 45 17 50

Philips Respironics United Kingdom
+44 800 1300 845

www.philips.com/respironics

For additional
information,
scan with your
smartphone.



Respironics, HomeLox, EverFlo and GoLox are trademarks of Respironics, Inc. and its affiliates.
All rights reserved.

Please visit www.philips.com/homelox



© 2011 Koninklijke Philips Electronics N.V.
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.
Geyer WMB 03/16/11 MCI 4104075 PN 1076174