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CONMED® 45° Angled ABC® Foot-control Pencil
 REF 130345
 Sterility guaranteed unless package has been opened, damaged or compromised.

Description:
 Angled (45°) Foot Control Disposable Handpiece with 10 foot (3.05m) cable.
 Allows the use of Monopolar ABC® Coagulation only.
 Can be used with ABC® Dissecting Electrodes.
 Sterile; Single Use

Intended Use:
 Handpiece is intended for use with the CONMED® ABC® Electrosurgical Generator for non-contact superficial ablation and hemostasis at the operative site.
 Maximum generator voltage 6.5V @ 3% duty cycle.

Contraindications:
 These devices should never be used when:
 ・ There is visible evidence of damage to the exterior of the device, such as cracks, cuts, punctures, nicks, abrasions, discoloration or connector damage.
 ・ Where conventional monopolar electrosurgery is inappropriate or unsafe, such as for small appendages, in circumcision or finger surgery.
 ・ In the presence of flammable gases, flammable prep solutions or drapes, oxidizing gases, such as Nitrous Oxide, or in oxygen-enriched environments.

Warning:
 ・ Venous gas emboli may be produced while utilizing Argon Beam Coagulation. To lower the risk of gas emboli, do not use Argon Beam Coagulation for prolonged periods, do not allow argon gas flow to make direct contact with open vessels, set the gas flow to the lowest setting capable of achieving the desired clinical effect, do not allow argon gas to flow for more than 5 seconds without cutting or coagulating tissue.
 ・ Do not directly secure the cable with metal instruments to surgical drapes. Activation of the device in contact with metal instruments may cause burns at the tissue/instrument interface.
 ・ Do not permit cables connected to associated electrosurgical accessories to be parallel to or in close proximity to the leads of other electrical devices.
 ・ To avoid burns never allow cable associated with this device to be in contact with skin of patient or touching operator.
 ・ To prevent burns, do not activate the electrosurgical unit if the tip is not in a position to deliver energy to target tissue.
 ・ Prolonged use of the handpiece causes the tip to become very hot and the tip remains hot for a period after use. Do not allow the tip to come in contact with the patient or others after activation to avoid burns.
 ・ **MODIFICATION OF THE DEVICE BY THE USER:** CONMED® warns explicitly against modifying the device, e.g., cutting off the tip of the malleable tube. Any modification voids warranty and exempts CONMED® from all liability.
 ・ Electrosurgery is potentially hazardous for patients with pacemakers, AICDs, neurostimulators or other active implants. High frequency energy may damage or adversely affect active implants. Ventricular fibrillation may occur.
 ・ Always place handpiece in a safe insulated location when not in use, to avoid burns.


Precautions:
 ・ Refer to CONMED® ABC® Generator operator manual and CONMED® ABC® Handpiece instructions prior to use.
 ・ **SINGLE PATIENT USE. DO NOT RESTERILIZE.** The ability to effectively clean and re-sterilize this device has not been established and subsequent re-use may adversely affect the performance, safety and/or sterility of the device.
 ・ Prevent contact of the argon tip to tissue.
 ・ Do not bend the shaft.
 ・ Use of the handpiece with insufficient argon gas flow may severely damage the tip.
 ・ Do not immerse or saturate handcontrol pencil or cord in fluid.
 ・ Do not pull on or stretch the handpiece cord.
 ・ Do not test the CONMED® by sparking the active electrode to the patient plate or other objects.
 ・ Use lowest possible power setting on the CONMED® ABC® Electrosurgical Generator capable of achieving the desired surgical effect.
 ・ Use shortest effective ABC® activation time necessary.
 ・ After patient use, this product should be considered a potential biohazard and disposed of per by hospital policy and applicable law.
 ・ The CONMED® ABC® Electrosurgical Generator should be used only by surgeons/healthcare practitioners experienced in Argon Beam electrosurgical techniques and safety.

Inspection
 These devices should be inspected before use. Visually examine the devices for obvious physical damage and do not use if found:
 ・ Cracked, broken or otherwise distorted plastic parts;
 ・ Broken or significantly bent handle, shaft or connector contacts.


Instructions for Use:

1. Before turning on generator power, insert the round connector end of the handcontrol pencil cord into the "Argon Beam Coagulator" receptacle of the ABC® Electrosurgical Generator and turn clockwise until tight.
 2. Plug three-pronged end connector of the compatible foot control switch into the "Beam Foot Control"-contacts of the ABC® electrosurgical generator.
 3. Remove tip protector prior to activation.
 4. Turn on generator power. To activate the pencil, depress footswitch. Argon gas will continue to flow approximately four seconds after footswitch is released.
 5. Turn off generator power before dismantling/ disposing of handpiece.


CAUTION: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.


© CONMED Corporation, 米国で出版


© CONMED Corporation - Gedrukt in de VS

© CONMED Corporation, Tryckt i USA.

© CONMED Corporation, 英国印制。

© CONMED Corporation, Printed in U.S.A.