Danish

ABC®Flex sonder

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English

ABC®Flex Probes

Sterility guaranteed unless package has been opened, damaged or compromised.

Description:
ABC®Flex Probes are designed for non contact superficial ablation and hemostasis through a flexible endoscope to the operative site. Probe Tip Indicator: A grey band, located 10 mm proximal from the tip of the ABC®Flex Probe, provides a visual reference point, from an endoscopic viewpoint, for probe tip positioning during the electrosurgical procedure. These probes should be used with the ConMed ABC® Electrosurgical Generator Systems 7500/7550 ONLY. These devices are compatible with endoscopes having an inner channel diameter of 2.5mm or larger. A compatible foot switch must be used to activate the device.
133023 = 7 Fr (2.3mm) x 220 cm
133270 = 7 Fr (2.3mm) x 270 cm
Sterile / Single Use

Intended Use:
The ABC® Flex Probes are sterile, single-use electrosurgical accessories used in conjunction with an ABC® generator for the delivery of argon gas and electrosurgical current through a flexible endoscope to perform coagulation.
Maximum output: 4.3kV, 80 W.

Contraindications:
These devices should never be used when:
• There is visible evidence of damage to the exterior of the device, such as cracks, cuts, punctures, nicks, abrasions, discoloration or connector damage.
• Where conventional monopolar electrosurgery is inappropriate or unsafe, such as for small appendages, in circumcision or finger surgery.
• In the presence of flammable gases, flammable prep solutions or drapes, oxidizing gases, such as nitrous oxide, or in oxygen-enriched environments.

Warnings:
POTENTIAL FOR COMBUSTION:
• **GI Tract:** Before using the ABC® device, ensure that no endogenous gases are present inside the bowel or around the surgical site; the bowel should be purged with nonflammable gas prior to ABC® electrosurgery.
• **Tracheobronchial Tree:** In the presence of oxygen or other flammable gases, material in the bronchoscope or tracheal tube can ignite during ABC® electrosurgery:
• Do not activate ABC® during O2 ventilation; O2 ventilation and ABC® activation must be alternated.
• Always keep the probe tip visible in the endoscope during activation to ensure proper position and performance.
• To avoid submucosal emphysema, do not activate the probe tip while in direct contact with the tissue.
• Avoid any contact of probe with **esophageal stents**.
• **Venous gas emboli may** be produced during ABC® electrosurgery.
• To lower the risk of gas emboli, do not use ABC® for prolonged periods; do not allow argon gas to flow for more than 5 seconds without coagulating tissue.
• Do not allow argon gas flow to make direct contact with open vessels.
• Set the gas flow to the lowest setting and shortest time capable of achieving the desired clinical effect.
• Monitor patient for early signs of pulmonary artery/venous emboli.
• Inadvertent burn avoidance:
• Do not activate the electrosurgical unit if the Flex Probe tip is not in a position to deliver energy to target tissue.
• Never allow the cable associated with this device to be in contact with skin of patient or touching operator. Do not permit cables connected to associated electrosurgical accessories to be parallel to or in close proximity to the leads of other electrical devices.
• The probe tip should extend **at least 10 mm** beyond the end of the endoscope; probe tip position can be monitored by observing the grey band (probe tip indicator) through the endoscope.
• Prolonged use of the probe causes the tip to become very hot and the tip remains hot for a period after use.
• Always place probe in a safe insulated location when not in use.
• Activate the foot switch while advancing the probe toward the tissue. Beam initiation will occur before contact with tissue. Do not allow probe tip to contact tissue.
• Ensure that the footswitch is not inadvertently depressed to prevent the probe from being unintentionally activated. Place footswitches in location that necessitates deliberate action in order to activate the device.
• Electrosurgery is potentially hazardous for patients with pacemakers, AICDs, neurostimulators and other active implants. High frequency energy may damage or adversely affect active implants. Ventricular fibrillation may occur.
• **MODIFICATION OF THE DEVICE BY THE USER:** CONMED® warns explicitly against modifying the device, e.g., cutting off the tip of the probe. Any modification voids warranty and exempts CONMED® from all liability.

Precautions:
• Refer to **CONMED® 7550/7500™ Electrosurgical Generator + ABC® Modes, Operator's Manual** instructions prior to use.
• **Single Patient Use - Do Not Resterilize:** The ability to effectively clean and re-sterilize this device has not been established and subsequent sterilization processing may adversely affect the performance, safety and/or sterility of the device.
• Use of the probe with insufficient argon gas flow may severely damage the tip.
• Do not immerse or saturate probe in fluid.
• Do not pull on, stretch, severely bend or kink the probe.
• After patient use, this product should be considered a potential biohazard and disposed of per hospital policy and applicable law.
• The ABC® Electrosurgical Generator should be used only by surgeons/healthcare practitioners experienced in Argon Beam electrosurgical techniques and safety.

Inspection:
These devices should be inspected before use. Visually examine the devices for obvious physical damage and replace if any defect is found:
• Broken, cracked, nicked, punctured, abraded, bent, distorted, discolored insulation / cable, connector or tip.

Instructions for Use:
1. Before turning on generator power, insert the large connector end of the probe into the "ABC" receptacle of the System 7500/7550™ electrosurgical generator.
2. Plug appropriate foot switch controller into generator.
3. Turn generator power on.
4. Verify that ABC is set to "Endo" mode. Set argon flow rate from 0.5 to 2 SLPM. **Do not exceed 2 SLPM.**
5. To activate the argon beam, operate the foot switch. Argon gas will continue to flow for approximately four seconds after button is released.
6. Turn off generator power before disconnecting accessories.

USA Federal law restricts this device to sale by or on the order of a physician.

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© CONMED Corporation, Trykt i USA.
© CONMED Corporation, Tryckt i USA
MDSS GmbH