## **Leather Care Instruction**

### A: Protect against excessive sunlight, dryness and heat.

- 1. Sunlight will fade and dry out leather.
- 2. Dryness will cause leather to crack.
- 3. Heat will dry out the leather and cause further damage.

### **B:** How to remove spots and stains:

- 1. Blot excess liquid up as quickly as possible with a clean 100% cotton white absorbent cloth or sponge.
- 2. If necessary, use clean lukewarm water only and let dry naturally. Do not dry with hair dyer, etc.
- 3. If water is used, clean the entire area where spot occurred seat cushion, arm, etc.

### Do not rub.

# C: How to remove butter, oil or grease:

- 1. Wipe excess butter, oil or grease off the leather with a clean 100% cotton white dry cloth. Then leave it alone as the spot should dissipate into the leather in a short period of time.
- 2. Do not apply water or try to wash the butter, oil or grease spot.

### **D:** How to remove chewing gum:

Place ice cubes in a plastic bag and apply on the chewing gum. Wait several minutes, lift the ice bag and remove the chewing gum.

#### E: General Care:

DO NOT USE CHEMICALS, SADDLE SOAP, CLEANING SOLVENTS, FURNITURE POLISH, OILS, VARNISH, ABRASIVE CLEANNERS, SOAP OR AMMONIA WATER.

Vacuum as a routine cleaning. Keep moist. Humidify if climate is dry.

These are recommended or suggested methods of cleaning. The manufacturer is not responsible for damage incurred while cleaning. Always try the cleaning method in a hidden area first to see results.