


REPLACEMENT PARTS



Product diagrams are for illustrative purposes only. Please confirm the correct parts needed prior to ordering by checking the product diagram online at: drivemedical.com or by contacting our Technical Service Team at techsupport@drivemedical.com

- | |
|---|
| 1- Back Only, Blue Ridge (D574BACK-BR), Jade (D574BACK-J), Rosewood (D574BACK-R) |
| 2- Head Rest Adjustment (D574P-1048) |
| 3- Push Handle (D574P-1050a) |
| 4- Lock Bar for Back (D574P-1050) |
| 5- Seat Frame (D574P-1053) |
| 6- Seat Only, Blue Ridge (D574SEAT-BR), Jade (D574SEAT-J), Rosewood (D574SEAT-R) |
| 7- Legrest, Blue Ridge (D574P-1055-BR), Jade (D574P-1055-J), Rosewood (D574P-1055-R) |
| 8- Footrest, Blue Ridge (D574P-1058-BR), Jade (D574P-1058-J), Rosewood (D574P-1058-R) |
| 9- Footrest Support Bar (D574P-1060) |
| 10- Recliner Mechanism, Left & Right (D574P-1035) |
| 11- Table Tray (D574P-1038) |
| 12- Front Caster (D574P-1027) |
| 13- Recliner Frame (D574P-1068) |
| 14- Rear Locking Caster (D574P-1025B) |
| 15- Armrest, Blue Ridge (D574P-1070-BR), Jade (D574P-1070-J), Rosewood (D574P-1070-R) |
| 16- Recliner Hardware Pack (D574P-1003) |
| 17- Side Panels, Left & Right (D574P-1020) |

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.