

250 LB
WEIGHT CAPACITY

3 Position Recliner

- D574-BR 3 Position Recliner, Blue Ridge, 1/cs
- D574-J 3 Position Recliner, Jade, 1/cs
- D574-R 3 Position Recliner, Rosewood, 1/cs

BARIATRIC
500 LB
WEIGHT CAPACITY

3 Position Recliner, Bariatric

- D574EW-BR 3 Position Extra Wide Recliner, Blue Ridge, 1/cs
- D574EW-J 3 Position Extra Wide Recliner, Jade, 1/cs
- D574EW-R 3 Position Extra Wide Recliner, Rosewood, 1/cs



D574-BR



AVAILABLE COLORS:



Rosewood Blue Ridge Jade

FEATURES:

- 3 positions: upright, deep recline and elevated legrest/footrest.
- Retractable lock bar secures the chair in the desired position.
- Self adjusting suspension in the seat prevents bottoming out.
- Side panels “pop-off” for easy cleaning and maintenance.
- Large, blow molded tray locks in 5 positions and can be stored on the side of the chair.
- 5” casters (2 with locks) make it easy to move and secure the recliner.
- Wide, comfortable armrests are angled to provide support in any reclined position.
- Neutral colored, sturdy steel frame is easy to clean and maintain.
- Moisture barrier on seat prevents seepage.
- Comfortable “built in” headrest.
- Foam and vinyl meet California Flammability. Standard and Boston Fire Department specifications.
- Limited Lifetime Warranty on Frame, 1 Year on Upholstery and Moving Parts.

SEE PARTS DIAGRAMS
ON PAGES 206-207

SPECIFICATIONS

Item #	Overall Height	Overall Width	Overall Depth	Width Between arms	Arm Height From Floor	Arm Height From Seat	Seat Height	Seat Width	Seat Depth	Carton Shipping Weight	Weight Limit
D574	47"	25"	41"	20"	26"	9"	21"	19"	19"	110 lbs.	250 lbs.
D574EW	47.5"	32"	41.5"	24.5"	26.5"	9.5"	21"	25"	20"	127 lbs.	500 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.