

# KANGA TS ADULT FOLDING TILT-IN-SPACE WHEELCHAIR

## PRICE LIST AND ORDER FORM



August 2014

Name: \_\_\_\_\_ Phone Number: \_\_\_\_\_

Account #: \_\_\_\_\_ PO #: \_\_\_\_\_ Date: \_\_\_\_\_

Bill To: \_\_\_\_\_ Marked For: \_\_\_\_\_

Address: \_\_\_\_\_ City: \_\_\_\_\_ Zip: \_\_\_\_\_

Ship to: \_\_\_\_\_

Address: \_\_\_\_\_ City: \_\_\_\_\_ Zip: \_\_\_\_\_

### KANGA TS – SEATING INCLUSIVE PRECONFIGURED\* AND READY TO SHIP VIA DRIVE MEDICAL RAPID DELIVERY PROGRAM - RDP

		MSRP
<input type="checkbox"/> <b>KG 1600</b>	Kanga TS, 16" Seating Included, Black Frame	\$3,700.00
<input type="checkbox"/> <b>KG 1800</b>	Kanga TS, 18" Seating Included, Black Frame	\$3,700.00
<input type="checkbox"/> <b>KG 2000</b>	Kanga TS, 20" Seating Included, Black Frame	\$3,700.00

\*Includes Planar Solid Seat, Planar Solid Backrest, Lateral Thoracic Supports, 90° Swing-away Footrests, Adjustable Angle Footplates, Heel/ Toe Loops, Headrest, Anti tippers, Removable Adductors & Abductor, Transit Tie Downs, Chest Harness, & Positioning Belt.

### KANGA TS – FRAME ONLY - INCLUDES TRANSIT TIE DOWN & ANTI-TIPPERS

		HCPCS	MSRP
<input type="checkbox"/> <b>KG A800</b>	Adult, Growable Frame, Black Frame Color*	E1161	\$2,550.00

\* Must order footrests separately.

### ACCESSORIES

SEAT PAN	HCPCS	MSRP	ELEVATING LEGREST	HCPCS	MSRP
<input type="checkbox"/> KG 1635 for use with KG A800 16"	E2231	\$245.00	<input type="checkbox"/> KG AELR Pair	E0990	\$235.00
<input type="checkbox"/> KG 1835 for use with KG A800 18"	E2231	\$245.00			
<input type="checkbox"/> KG 2035 for use with KG A800 20"	E2231	\$245.00			
<b>SOLID PLANAR SEAT</b>			<b>LATERAL THORACIC SUPPORTS</b>		
<input type="checkbox"/> KG 1633 Solid Seat, 16"	E2292	\$149.00	<input type="checkbox"/> KG 1039 Pair, for KG A800	E0956	\$159.00
<input type="checkbox"/> KG 1833 Solid Seat, 18"	E2292	\$149.00			
<input type="checkbox"/> KG 2033 Solid Seat, 20"	E2292	\$149.00			
<b>SOLID PLANAR BACK</b>			<b>HEEL AND TOE LOOPS</b>		
<input type="checkbox"/> KG 1632 Solid Back, 16"	E2291	\$195.00	<input type="checkbox"/> KG 8046 For KG A800	E0951/ E0952	\$149.00
<input type="checkbox"/> KG 1832 Solid Back, 18"	E2291	\$195.00			
<input type="checkbox"/> KG 2032 Solid Back, 20"	E2291	\$195.00			
<b>HEADREST</b>			<b>ADDUCTOR</b>		
<input type="checkbox"/> KG 8000 Headrest	E0955	\$198.00	<input type="checkbox"/> KG 1040 Removable, Pair	E0956/ E1028	\$120.00
<b>90° FOOTRESTS HANGERS WITH ADJUSTABLE ANGLE FOOTPLATES</b>			<b>ABDUCTOR</b>		
<input type="checkbox"/> KG A045 Pair, use with 16"/18"/20"	K0052/K0040	\$179.00	<input type="checkbox"/> KG 1041 Removable	E0957/ E1028	\$126.00
<b>ANTI TIPPERS</b>			<b>UPPER EXTREMITY SUPPORT</b>		
<input type="checkbox"/> KG 1020 Pair	E0971	\$62.00	<input type="checkbox"/> KG 8025 For use with KG A800	E0950	\$225.00
<b>ELEVATING LEGREST</b>			<b>CANOPY</b>		
<input type="checkbox"/> KG AELR Pair	E0990	\$235.00	<input type="checkbox"/> KG 8026 Use with KG P800, KG 1000 KG 1200, KG 1400	K0108	\$230.00
			<b>POSITIONING BELT</b>		
			<input type="checkbox"/> KG 2042 For use with KG A800	E0978/E0960	\$70.00

\* Included with frame only.

SPECIFICATIONS	KG A800	SPECIFICATIONS	KG A800
Overall Width (Open)	24"-28"	Back Cane Height	19"
Overall Width (Closed)	12"	Push Handle Adjustment Range (Floor to Handle)	39" - 44"
Seat Depth (Adjustable)	17"-19"	Seat Rail to Footplate Range	8" - 19"
Seat Width (Adjustable)	16"-18"-20"	Arm Height Range	10" - 13"
Product Weight	29 lbs.	Overall Length without Footrests	33"
Weight Capacity	300 lbs.	Overall Length with Footrests	42"
Seat to Floor Height	19"	Back Angle Adjustment Range	85° - 115°

Drive Medical Design and Manufacturing

99 Seaview Boulevard, Port Washington, NY | t: 877.224.0946 | f: 516.998.4601 | www.drivemedical.com