


SPITFIRE EX 1320-3 WHEEL

PARTS DIAGRAM



Product diagrams are for illustrative purposes only. Please confirm the correct parts needed prior to ordering by checking the product diagram online at: drivemedical.com or by contacting our Technical Service Team at techsupport@drivemedical.com

SPECIFICATIONS

ITEM #	SPITFIRE132016FS12, SPITFIRE132016FS21
HCPCS	GROUP 1 / K0800
WEIGHT CAPACITY	300 LBS.
TOP SPEED	4 MPH
TOP RANGE	9 MILES (12AH)/15 MILES (21AH)
TURNING RADIUS	46"
CLIMBING ANGLE	6°
GROUND CLEARANCE	3"
DIMENSIONS (L X W X H)	41" X 19" X 36"
FLOOR TO SEAT HEIGHT	20"- 23"
SEAT DIMENSIONS (W X D X H)	17" X 16" X 14"
CONTROLLER	PG 545
MOTOR/DRIVE SYSTEM	24V X 270W X 4700RPM
BATTERIES	12V X 12AH/21AH (2)
CHARGER	2 AMPS OFFBOARD
BRAKES	ELECTROMAGNETIC
FRONT WHEEL	7.5" X 2" FLAT FREE
REAR WHEELS	8.5" X 3" FLAT FREE
BASE WEIGHT	66 LBS.
BATTERY WEIGHT	20 LBS. (12AH) 28 LBS. (21AH)
SEAT WEIGHT	17 LBS.
TOTAL WEIGHT	111 LBS. (12AH) /121 LBS. (21AH)
HEAVIEST PIECE WEIGHT	34 LBS. FRONT SECTION
WARRANTY	LIFETIME ON FRAME, 14 MONTHS ON ELECTRONIC CONTROLLER AND DRIVE TRAIN COMPONENTS, SIX MONTHS ON BATTERY