

**Harmar®****SUMMIT STAIR LIFT**☐ **QUOTE** ☐ **ORDER**

Date_____

Account Number_____

Sales Representative_____

Purchase Order Number_____

Dealer_____

Ship To Address_____

Address_____

City_____ State_____ Zip_____

Phone_____

Fax_____

Contact_____


Email_____

Phone_____

PARTS & OPTIONS

Model Number	Rail Length	Custom Cut Rail
<input type="checkbox"/> SL350AC	<input type="checkbox"/> Standard Rail Length - 16'	<input type="checkbox"/> SLR16C Less than 16'
<input type="checkbox"/> SL350OD	<input type="checkbox"/> SLRU20 20' Rail Upgrade	<input type="checkbox"/> SLR20C 16' to 20'
<input type="checkbox"/> SL500HD	<input type="checkbox"/> SLRU25 25' Rail Upgrade* (250 lb capacity)	<input type="checkbox"/> SLR25C 20' to 25'* <i>Custom Cut part numbers include track. *If track length is greater than 20', maximum lift capacity is reduced to 250 lbs</i>

Upholstery	Rail Length	
<input type="checkbox"/> Almond	<input type="checkbox"/> SL310 Keylock on Lift Only	<i>*Not available on SL500HD unit. SL500HD unit includes 23" seat and pendant control.</i>
<input type="checkbox"/> Evergreen	<input type="checkbox"/> SL315 Key Lock Package with Call/Send	
<input type="checkbox"/> Mocha	<input type="checkbox"/> SL330 Grocery Basket* (mounts to side of seat) *N/A with SL325	
	<input type="checkbox"/> SL320 19" Folding Arms with Pendant*	
	<input type="checkbox"/> SL325 21" Wide Seat*	

HOW TO MEASURE**Measurements**_____ **A to B** - Inches from top landing to bottom floor_____ **A to H** - Inches from top landing to nearest obstruction—for example, a door**Rise** _____ inches **Tread** _____ inches☐ **Left** ☐ **Right****Notes**
