




342 Steel Frame Screen Parts List



PARTS LIST				
Image Number	Note	Part Number	Description	QTY per Product
1		240100	342 COMPLETE FRAME	1
2	A	211650	SCREEN LEG W/CASTERS ASSEMBLY KIT (313-3	2
2A	A	211610	DOUBLE LEG 313/342	2
2A	A	702220	CASTER INSERT, 1" RND X 16 - 18 GA	4
2A	A	400000	CASTER 2" BLACK/TWIN HOODED AP7	4
2B	A	566506	SCREW-MACH 1/4-20 X 1 1/2 RND PHIL M/S ZINC	4
2B	A	566204	NUT-1/4-20 HEX NYLON LOCKNUT ZINC	4
2E		701700	CRUTCH TIP 1" BLACK	4
NOTES				
Note	Description			
A	It is NOT recommended that these parts be sold seperately however in special cases they can be.			

342 Steel Frame Screen Parts List



PARTS LIST				
Image Number	Note	Part Number	Description	QTY per Product
1A		299700	VINYL PANEL 342 WHITE	1
1A	A	2997XX	342 SURE CHEK PANEL	1
1B		240115	SCREEN ROD 342	2
1C		704500	ROD & ROD SPRING RETAINER(SET)	4
NOTES				
Note	Description			
A	XX- For SURECHEK use color code for specific part number i.e. 93 - Mint (299793)			