





851 Treatment Table Parts List




PARTS LIST

Image Number	Note	Part Number	Description	QTY per Product
1	A	U3004XX	851 UPHOLSTERED TOP	1
1	B	222700	TABLE FRAME 850/851/855	2
1	B	701010	1X2X18 RECT INSERT CAP/GLIDE LT.GRAY	4
1	B	566519	SCREW-MACH 1/4-20 X 2 1/2 PAN PHILLIP M/S ZINC	4
1	B	566528	SCREW-TAP 10 X 2 3/4 PAN PHIL T/S ZINC	4
2		219100K	FACE PAPER HLDR ASSM. KIT	1
2	B	219100	FACE PAPER HLDR COMPLETE	1
2	B	566540	SCREW-TAP 6 X 1 BUGLE PHILLIP DRYWALL SCREW	2
3		219120K	FACE SLOT TUBE ASSM. KIT	1
3A	B	219120	FACE SLOT TUBE	1
3B	B	566526	SCREW-TAP 10 X 1 3/4 PAN PHIL T/S ZINC	2
3B	B	566315	WASHER-1/4 NYLON SPACER (.257 X .5 X 1.0)	2
4		219165K	851 FACE PAPER CUTTER ASSM. KIT	1
4	B	219160	FACE PAPER CUTTER	1
4	B	705600	DOUBLE SIDED FOAM TAPE, 3/4" X 36YDS	8 in.

NOTES

Note	Description
A	XX- Use color code for specific part number i.e. 02 - Royal Blue
B	These parts are included with the UPHOLSTERED TOP

851 Treatment Table Parts List



PARTS LIST				
Image Number	Note	Part Number	Description	QTY per Product
1		222708K	WELDED LEG ASSM. KIT 19"	2
1	A	222708	WELDED LEG ASSM. 19"	1
1		222709K	WELDED LEG ASSM. KIT 25" (850 style)	2
1	A	222709	WELDED LEG ASSM. 25" (850 style)	1
1		222710K	WELDED LEG ASSM. KIT 30" (850 style)	2
1	A	222710	WELDED LEG ASSM. 30" (850 style)	1
1		222711K	WELDED LEG ASSM. KIT 36" (850 style)	2
1	A	222711	WELDED LEG ASSM. 36" (850 style)	1
1	A	500350	GLIDE NON-SWIVEL ZINC 5/16X1" NP LN	2
1	A	500105	INSERT, BLK HDPE, RIBBED 16-18 ga 1X2 , 5/16-18 THREAD	2
1A	A	566506	SCREW-MACH 1/4-20 X 1 1/2 RND PHIL M/S ZINC	2
1B	A	566519	SCREW-MACH 1/4-20 X 2 1/2 PAN PHILLIP M/S ZINC	2
1B	A	566204	NUT-1/4-20 HEX NYLON LOCKNUT ZINC	2
2A		222705K	PAPER ROD HOLDER ASSM. KIT (850 STYLE)	1
2A	A	222705	850/851/855/856 ROD, PAPER HOLDER 1/2" ALUM	1
2A	A	700650	1/2" PLASTIC CAP, WHITE	2
2B		201810	PAPER ROD CLAMP, GRAY	2
2C		271500	VINYL TABLE STRAP	1
2D		566302	WASHER-#10 SAE F/W ZINC	2
2D		566534	SCREW-TAP 10 X 3/4 RND WASHER SQ T/S ZINC	2
3		222720	850/851/855/856 WELDED BRACE ASSM.	1
3		566401	BOLT-CARRIAGE 1/4-20 X 1 1/2 ZINC	2
3		566204	NUT-1/4-20 HEX NYLON LOCKNUT ZINC	2
4	B	222730	WELDED, 850-851-856 REINF. BRACE ASSM WELDED	1
4		566401	BOLT-CARRIAGE 1/4-20 X 1 1/2 ZINC	2
4		566204	NUT-1/4-20 HEX NYLON LOCKNUT ZINC	2
NOTES				
Note	Description			
A	These parts are included in ALL above listed kit.			
B	For tables with optional shelf, cabinet, or drawer. Must order double-sided tape listed on Page 3			

851 Treatment Table Parts List



PARTS LIST				
Image Number	Note	Part Number	Description	QTY per Product
1		263960K	850 SERIES PLAIN SHELF KIT	1
1	A	263960	850 SERIES PLAIN SHELF	1
1	A	705600	DOUBLE SIDED FOAM TAPE, 3/4" X 36YDS	1 ft.
2		263950	850,851,856 MODULE SHELF	1
3		266780K	850/851/856 CABINET MODULE ASSM. KIT	1
3A	A	266780	850/851/856 CABINET MODULE	1
3A	A	566519	SCREW-MACH 1/4-20 X 2 1/2 PAN PHILLIP M/S ZINC	4
3A	A	566204	NUT-1/4-20 HEX NYLON LOCKNUT ZINC	4
4		266760K	850/851/856 MODULE DRAWER ASSM. KIT	1
4	A	266760	850/851/856 MODULE DRAWER ASSM	1
4A	A	566519	SCREW-MACH 1/4-20 X 2 1/2 PAN PHILLIP M/S ZINC	4
4A	A	566204	NUT-1/4-20 HEX NYLON LOCKNUT ZINC	4
NOTES				
Note	Description			
A	These parts are included in the above listed kit.			