

Wrist Type Automatic  
Blood Pressure Monitor



## Main Components

Model: LD732



Acryl Faceplate

LCD Display

Battery Cover

O/I Button

Memory Button

Wrist Cuff

## Key Features

Model: LD11A

## Size & Weight



Weight: 110g without batteries

Dimension: 70mm(L) x 60mm(W) x 30mm(H)

## Key Features

**946-01-BLU**

One touch operation to finish the measurement



The battery cover is easy to be opened by pressing the two sides

Concave design





## Key Features

Model: LD732

## Operation



Standby Mode



LCD Full Display



'0' will appear on  
the screen



Automatic Inflation



Automatic Deflation



The results with measuring time display  
together in two screens alternately


When the measurement is finished, the value is recorded automatically

## WHO Indication

### Key Features

**946-01-BLU**

### WHO Blood Pressure Indication



Grade	Sys.	Dia.
Optimal	< 120mmHg	< 80mmHg
Normal	120 ~ 129mmHg	80 ~ 84mmHg
High-normal	130 ~ 139mmHg	85 ~ 89mmHg
Mild Hypertension	140 ~ 159mmHg	90 ~ 99mmHg
Moderate Hypertension	160 ~ 179mmHg	100 ~ 109mmHg
Severe Hypertension	180mmHg	110mmHg

Note: If the systolic blood pressure and diastolic blood pressure fall into different categories, the higher value should be taken for classification.

## Key Features

**946-01-BLU**

The monitor provides a blood pressure and pulse rate measurement even when an irregular heartbeat occurs.

**Note 1:** It is important that you be relaxed, remain still and do not talk during measurement.

**Note 2:** We recommend contacting your physician if you see this “  
indicator frequently.



## IHB Detector



## Detailed Specs

**946-01-BLU**

Measuring Method	Oscillometry
Measuring Range	40 to 260 mmHg (blood pressure)
	40 to 160 beats/minute (pulse rate)
Measuring Accuracy	$\pm$ 3 mmHg for static pressure
	$\pm$ 5% of the reading for the pulse rate
Inflation	Automatic
Rapid Deflation	Electromagnetic valve
Power Source	2 “AAA” $\times$ 1.5V
Memory	90 sets of memories
Operation Temperature and Humidity	+ 10°C to + 40°C, 85% and below
Storage Temperature and Humidity	-20°C to + 50°C, 85% and below
Wrist Size	Applicable for wrist size 12.5 cm-20.5 cm
Monitor Size	70mm(L)*60mm(W)*30mm(H)
Weight	Approximately 110g without batteries