

Afiscooter-C3

The outstanding design of the Afiscooter-C along with the user comfort facilities:

A fully adjustable swiveling orthopedic seat that gives easy access and an enjoyable ride, the DELTA tiller that provides easy steering and is fully adjustable to accommodate any driver. A top quality drive is guaranteed with an advanced full suspension system. In addition, the Afiscooter-C features a large rear basket. The 3 wheel model allows for perfect mobility on narrow sidewalks and convenient maneuvering in shopping centers for both indoor and outdoor use.



Afiscooter-C3



Afiscooter-C3
Metallic Blue



Afiscooter-C3
Metallic Silver

Colors available for these models 

Afiscooters

Afiscooters develops and manufactures advanced, electrically-powered mobility scooters since 1978. Afiscooters' mobility scooters are designed with reliability and driver safety in mind. Our worldwide distribution network offers year-round local technical support. Afiscooters' scooters meet European ISO 9001-2000 and EN-12184 certification and US safety standards with FDA approval.

1-800-809-3010 (US only)

www.afiscooters.com

Features

- ▶ A fully adjustable orthopedic seat provides comfort while driving and swivels for ease of boarding
- ▶ Easy tiller adjustment for any driver
- ▶ Fits standard car lifts for the traveling user
- ▶ An optional canopy is available for sun protection



All operational controls including lights, horn and battery display are conveniently located in the center of the one-hand operation Delta handlebar



Plenty of storage space in the standard: front basket and optional large rear basket



Specifications

		Afiscooter-C3	
Overall length		51"	130 cm
Overall width		26"	67 cm
Overall height	Without canopy	50"	126 cm
Ground clearance		4"	10 cm
Weight	With batteries	225 lb	102 kg
	Without batteries	161 lb	73 kg
Travel range		25 m	40 km
Weight capacity Inc. driver		330 lb*	135 kg
Maximum speed		9.3 mph	15 Km/h
Maximun climbing angle		16 %	
Turning radius		50"	128 cm
Seating		18"	46 cm
Tires type		Pneumatic	
Tires size	Front	4.1/3.5X6"	
	Rear	4.1/3.5X6"	
Battery type		Two 12V 45-50A/h	
Power unit		Rear wheel drive, 24V DC motor	
Color		Metallic Silver / Blue / Cherry Red	

* US Model



AFIKIM
ELECTRIC VEHICLES

1-800-809-3010
www.afiscooters.com