

# Caddy

The world's first folding scooter. The Caddy allows for compact storage at home, transport in the trunk of the family car or on a plane. Featuring a strong light weight aluminum chassis and easy front-wheel drive, the vehicle provides maximum maneuverability and durability. It takes just a few seconds to fold or open the Caddy which makes it the ideal scooter for trips abroad.



Quick folding for transport in the trunk of a car

Scan the code to watch the video >>



## AfiScooters

AfiScooters develops and manufactures advanced, electrically-powered mobility scooters.

AfiScooters' mobility scooters are designed with reliability and driver safety in mind. Our worldwide distribution network offers year-round local technical support. AfiScooters' scooters meet European ISO 9001-2000 and EN-12184 certification and US safety standards with FDA approval.

1-800-809-3010 (US only)

[www.afiscooters.com](http://www.afiscooters.com)

## Features

Easily disassembled into five light-weight parts



Easy front-wheel drive



Easy and simple operation



## Specifications

		Caddy	
Overall length		48"	121 cm
Overall width		24"	61 cm
Overall height		38"	97 cm
Ground clearance		3"	8 cm
Weight	With batteries	93 lb	42 kg
	Without batteries	62 lb	28 kg
	Weight of heaviest part	37 lb	17 kg
Travel range		12 m	20 km
Weight capacity Inc. driver		220 lb	100 kg
Maximum speed		4 mph	6 Km/h
Maximun climbing angle		6%	
Turning radius		41"	105 cm
Seating		19"	48 cm
Tires type		Full Tire	
Tires size	Front	250X50 mm	
	Rear	250X50 mm	
Battery type		Two 12V 17A/h	
Power unit		Front wheel drive, 24V DC motor	
Color		Blue	



**AFIKIM**  
ELECTRIC VEHICLES

1-800-809-3010  
[www.afiscooters.com](http://www.afiscooters.com)