

Afiscooter-M

Afiscooter-M is a powerful electric vehicle designed for service and maintenance personnel due to its impressive carrying capacity and maneuvering capabilities.

Afiscooter-M is a super strong vehicle with an advanced ergonomic design to handle loads up to 300 kg.

Afiscooter-M is equipped with an extremely large loading surface, long travel range and superb traction. This vehicle allows particularly easy steering and a smooth driving experience, in a variety of environments: corridors, service elevators, warehouses and narrow passages. The independent suspension system and the orthopedic chair guarantee that you will be comfortable during a full working day.



Afiscooter-M



Afiscooter-M
Extended loading surface



Afiscooter-M
Loading surface with canopy

Afikim Electric Vehicles develops and manufactures advanced, electrically-powered mobility scooters since 1978. Afiscooters' mobility scooters are designed with reliability and driver safety in mind. Our worldwide distribution network offers year-round local technical support. Afikim Electric Vehicles meet European ISO 9001-2000, EN-12184 certification, and US safety standards with FDA approval.

1-800-809-3010 (US Only)

www.afiscooters.com


Specifications

Afiscooter-M

Overall length		89"	225 cm
Overall width		35"	90 cm
Overall height	Without canopy	49"	125 cm
	With canopy	65"	164 cm
Ground clearance		5"	12.5 cm
Rear Wheel Track	Inner	20"	50 cm
	Outer	35"	88 cm
Wheel Base		76"	192 cm
Weight	With batteries	463 lb	210 kg
	Without batteries	300 lb	136 kg
	Add for canopy	+26 lb	+12 kg
	Add for large platform		+3 kg
Travel range		33 m	50 km
Rear platform length	Short/Long	37" / 45"	95 cm / 115 cm
Rear platform width		35"	88 cm
Weight capacity Inc. driver		661 lb	300 kg
Maximum speed		7.5 mph	12 kmh
Maximum climbing angle		13%	
Turning radius		79"	200 cm
Seating		18"	46 cm
Tires type		Pneumatic	
Tires size	Front	3X10"	
	Rear	16X7.5-8.0"	
Battery type		Two 12V 150Ah	
Power unit		Rear wheel drive, 24V DC motor	
		1400w 2500w MAX	
Controler type		SD 200 A	
Charger type		8A	
Color		Metallic Silver	

Advantages

- Bright, easy to read LCD display control panel indicates speed, driving length, working hours' counter, battery status, lighting and signing.



- Designed sun/rain canopy as an integral part of the vehicle, for working smoothly in all weather conditions. Rain cover is optional.



- Optional: Standard loading surface / Extended loading surface / 25 cm walls and a drop-down tailgate.



AFIKIM
ELECTRIC VEHICLES

1-800-809-3010
www.afiscooters.com