



Toxicology

Dependable results and simple to administer.

# *i*Cassette<sup>®</sup>

Urine Drug Screening Device.




Substance abuse testing with more **substance**.

## How it works.

PRODUCT TRAINING AND  
CERTIFICATION AVAILABLE ONLINE.

### Features & Benefits

- » Rapidly screen up to 11 illicit and prescription drugs
- » FDA 510(k) Cleared; CLIA Waived 6 drug test available for clinical use
- » Simple procedure: collect and pipette to test—results in as little as 5 minutes
- » Test device remains in control of administrator
- » Long shelf-life



PRODUCT PROCEDURE NOTE: Refer to product insert for complete instructions, limitations, and warnings.

For in vitro diagnostic use only. This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgement should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

### Urine Screening Device Cutoff Levels

Test	Calibrator	Cutoff
Amphetamine (AMP)	d-Amphetamine	1,000 ng/mL
Barbiturates (BAR)	Secobarbital	300 ng/mL
Benzodiazepines (BZO)	Oxazepam	300 ng/mL
Buprenorphine (BUP)	Buprenorphine	10 ng/mL
Cocaine (COC)	Benzoylcegonine	300 ng/mL
Marijuana (THC)	11-nor- $\Delta^9$ -THC-9-COOH	50 ng/mL
Methadone (MTD)	Methadone	300 ng/mL
Methamphetamine (mAMP)	d-Methamphetamine	1,000 ng/mL
Methylenedioxymethamphetamine (MDMA) Ecstasy	d,l-Methylenedioxymethamphetamine	500 ng/mL
Opiate (MOP)	Morphine	300 ng/mL
Opiate (OPI)	Morphine	2,000 ng/mL
Oxycodone (OXY)	Oxycodone	100 ng/mL
Phencyclidine (PCP)	Phencyclidine	25 ng/mL
Propoxyphene (PPX)	Propoxyphene	300 ng/mL
Tricyclic Antidepressants (TCA)	Nortriptyline	1,000 ng/mL

