

ASPIRATOR

Model 280



MODEL 280

ASPIRATOR – The Gomco® Model 280 Aspirator is designed for general suction use in the home, hospitals or clinics. It is compact, reliable and economical.

FEATURES

- Ideal for homecare use
- Enclosed cabinet protects motor chamber
- Low maintenance diaphragm pump
- Hospital grade electric plug
- Lightweight for easy portability
- Quiet operation
- Backed by Gomco's reputation for quality and performance
- Can be ordered with plastic disposable collection unit or glass bottle

MODEL 280 – The heart of the Model 280 is Gomco's diaphragm pump, designed to provide low maintenance and longlasting, quiet performance. The pump is designed for greater dependability, requires no oiling and is completely enclosed for maximum protection.

The Model 280 incorporates Gomco's High Efficiency Bacteria Filter, a hydrophobic, microporous membrane which filters air with maximum efficiency while helping to block the flow of aqueous fluids and aerosol contaminants. The unit has a lighted ON/OFF switch, hospital grade electrical plug for added user protection, a 600 ml collection bottle or an 1100 ml disposable canister, and a 6-foot length of disposable tubing. Model 280 is UL listed.

A built-in handle and its light weight allows transporting the Model 280 with minimal effort.

MODEL	COLLECTION BOTTLE	MOTOR	COMPRESSOR	VACUUM RANGE	FLOW RATE	DIMENSIONS HxDxW	SHIPPING WEIGHT
280	Glass	1/16 hp	Diaphragm	0-22 in. Hg	12 lpm	10.5 x 8 x 11.5 in. 26.7 x 20.3 x 29.2 cm	18 lbs. 8 kg
2801	Disposable	1/16 hp	Diaphragm	0-22 in. Hg	12 lpm	10.5 x 8 x 11.5 in. 26.7 x 20.3 x 29.2 cm	18 lbs. 8 kg

ORDERING INFORMATION

w/Glass	w/Disposable	Voltage
01-22-0280	01-22-2801	115/60
01-23-0280	01-23-2801	220/50
01-24-0280	01-24-2801	220/60