

# Disposable Collection Canisters



## Discover the Superior Performance of Allied's Enhanced Design.

Superior design and engineering make the Allied Disposable Collection Canister more efficient to work with than any others you may have tried ... making the job of replacing and disposing of collection canisters a lot easier.

- The snap-on lid is easy to attach to the canister and virtually impossible to remove ... and the audible “snap” gives clear confirmation that the lid is properly sealed and secure.
- With our unique DISS connection accessory, Allied canisters can be attached directly to most portable floor suction units and all wall suction units in the OR, ER, ICU and CCU departments. This helps eliminate excess tubing and potentially troublesome wall brackets.
- A filter located in the lid is designed to enhance performance over previous float systems. The filter design and self-sealing material improve flow 20%. When the filter comes into contact with moisture, it immediately stops flow from the patient and prevents fluids from entering the system.
- Built-in tubing hooks in the lid eliminate the need to wrap excess tubing around the regulator, preventing kinks and occlusions. Caps for all ports are attached to the lid for convenience in sealing the canister after use.
- The canister label provides accurate, easy-to-read calibrations and an enlarged area in which to write patient information.





*The unique DISS connection attaches directly to most suction units, saving time and eliminating wall brackets*



#### Canister Holder

Holds collection canister by inserting into slide bracket.  
20-02-0150



#### Collection Canister Carrier

Carries up to four disposable collection canisters of either 1500 ml or 2400 ml size. Easy to clean, durable chrome-plated steel construction. No-snap, free-spinning rotating casters. Comes with a DISS swivel-nut connection for attaching a suction regulator, a 5/16" barb fitting for attaching a vacuum hose, and two shut-off valves for regulating the suction to the canister. Four slide brackets are permanently affixed to the pole. 20-02-0023

#### Slide Bracket Kit

Attaches to 20-02-0031 collection canister carrier pole to accommodate four additional canister holders. Kit includes four canister holders, four slide brackets, and self-tapping screws. 20-02-0028



#### Specimen Bag

Filters large particles from aspirate. Fits into accessory port for easy accessibility. 4.5". 20-02-0151



#### Critical Measurement Insert

For precise volume determination of aspirated material. Inserts into 1500 ml canister. 20-02-0152

### ORDERING INFORMATION

Disposable Containers	Size	Inlet	Units/Case	Accessories	Units/Carton	
20-08-0001	2400 ml	DISS	36	20-02-0023	Collection Canister Carrier,	
20-08-0002	2400 ml	Stem	36		4-canister model, 36" H	1
20-08-0003	1500 ml	Stem	48	20-02-0030	Collection Canister Carrier,	
20-08-0004	1500 ml	DISS	48		4-canister model, 21" H	1
20-08-0005	2400 ml	DISS w/tubing*	36	20-02-0031	Collection Canister Carrier,	
20-08-0006	2400 ml	Stem w/tubing**	36		8-canister model, 52" H	1
20-08-0007	1500 ml	Stem w/tubing**	42	20-02-0028	Slide Bracket Kit	4
20-08-0008	1500 ml	DISS w/tubing*	42	20-02-0150	Canister Holder	12
*PVC Tubing, 1/4" ID, 72"				20-02-0151	Specimen Bag	12
**PVC Tubing, 1/4" ID, 18" & 72"				20-02-0152	Critical Measurement Insert	12

\*PVC Tubing, 1/4" ID, 72"

\*\*PVC Tubing, 1/4" ID, 18" & 72"



**ISO 9001CERTIFIED**

Allied Healthcare Products, Inc. • 1720 Sublette Avenue • St. Louis, Missouri 63110  
314/771-2400 • 800/444-3954 • Fax: 314/771-0650 • [www.alliedhpi.com](http://www.alliedhpi.com)  
International offices in Europe, Asia, and South America

U.S.

Toll Free: 800/444-3940  
Fax: 800/477-7701

Canada & Mexico

Toll Free: 800/446-0552  
Fax: 800/246-6201

International

Telephone: 314/268-1683  
Fax: 314/771-5183

©Allied Healthcare Products, Inc. All specifications are nominal and subject to change without notice Form 22-00-0004 Rev. 5/04