

Dispette2



The standard in ESR testing.

Economic
Precise
Safe
Reliable

The Dispette2 is a CLIA waived test that meets the Westergren Standard. It is available in a convenient 100 count box or in bulk to meet your needs.



Dispette2



- Precision manufacturing of the pipette ensures a constant bore size which gives improved reproducibility of results.
- Unique "triple plug" tip design provides enhanced protection barrier from accidental over-filling of the pipette and the spilling of blood.
- Pierceable reservoir cap provides greater protection against blood spillage.
- "Sealing" ring molded into the reservoir below the fill line ensures easy filling of Dispette tube.
- Restricted end provides a buffered stop to the downward progress of the Dispette indicating that it is near to its final position.
- Reservoirs are pre-filled with 0.25 mL of 0.9% saline ensuring the correct dilution to comply with the standard of Westergren ESR testing when 1 mL of EDTA blood is added.
- Enlarged chamber between the reservoir cap and base provides an improved mixing capacity for the blood sample and the pre-filled diluent.

WESTERGREN

Westergren pipette dimensions together with the correct dilution of the blood sample are the internationally accepted standard for ESR testing.

The Westergren "stamp" on the Dispette 2 authenticates the product as conforming to this standard, which is recommended by both the Clinical Laboratory Standards Institute in the USA and the International Council for Standardization in Haematology.

Dispette is a registered trademark of Guest Scientific.

