



**AM-680 Training Stairs with Bus Step**

Maple plywood base and handrails with natural finish. Three 4" risers on one side and two 6" risers on the other. When the bus step is installed, a 12" riser is created. The platform is 24" x 30".

Handrail Height 36"; Overall Height 48"  
450 lb. Capacity; Footprint: 55" x 31"

**AM-686 Mobile Posture Mirror**

Maple plywood with natural finish  
23" x 70" with 3" heavy duty Caster

Height: 76"  
Width: 26"  
Depth: 18"



**AM-840 Footstool**

All welded steel construction with silver vein powder coating. Stepping surface / platform 10" x 15" with non-skid rubber tread.

Overall size 14"W x 18"L x 10"H  
300 lb. Capacity



**AM-842 Footstool with Handrail**

All welded steel construction with silver vein powder coating. Stepping surface / platform 10" x 15" with non-skid rubber tread.

Overall size 14"W x 18"L x 10"H  
Step with handrail height of 41".  
300 lb. Capacity



**AM-844 Mobile Stool**

All welded steel construction with silver vein powder coating. Seat is 14" x 14" with 2" firm density foam. Height is 18" with 4 swivel casters.

225 lb. Capacity; 2.5" casters

**AM-846 Pneumatic Stool**

Seat is 15" diameter with 3" of Firm Density foam padding. Seat height is adjustable from 17" to 22". Base has 5 legs and is 21" overall width.

300 lb. Capacity



**Other Options and Accessories:**

- 4 Caster Base
- 6" Clearance Base
- Face Plug
- AM-805 Paper Holder with Cutter
- AM-807 Fold Down Side Rails
- Face Cut Out