



AM-BA250 Treatment Table

- | | | | |
|--|---|--|---|
| <ul style="list-style-type: none"> • 2-Section Treatment Table • 3 Piece Head • 27"W x 76"L | <ul style="list-style-type: none"> • Min. Height: 18" • Max Height: 37" • 1 1/2" Firm Density Foam | <ul style="list-style-type: none"> • Head: 15" x 13" with Contoured Face Opening • Arms: 12.5" x 5.5" • Body: 61" x 27" | <ul style="list-style-type: none"> • 400 lb. Capacity • 4 Total Lock Swivel Casters |
|--|---|--|---|



AM-BA350 Treatment Table

- | | | | |
|--|---|---|---|
| <ul style="list-style-type: none"> • 3-Section Treatment Table • 27"W x 76"L | <ul style="list-style-type: none"> • Min. Height: 18" • Max Height: 37" • 1 1/2" Firm Density Foam | <ul style="list-style-type: none"> • Head: 15" x 27" with Contoured Face Opening • Center: 22" x 27" • Foot: 39" x 27" | <ul style="list-style-type: none"> • 400 lb. Capacity • 4 Total Lock Swivel Casters |
|--|---|---|---|