

# AM-BAX SERIES

BAR ACTIVATED HI-LO CONTROL



## AM-BAX5000 Treatment Table with Motorized Center Section

- 3-Section Treatment Table
- 3 Piece Head
- 27"W x 72"L
- Min. Height: 19"
- Max Height: 38"
- 1 1/2" Firm Density Foam
- Head: 15" x 13" with Contoured Face Opening
- Arms: 15" x 7"
- Center: 19" x 27"
- Foot: 38" x 27"
- 400 lb. Capacity
- 4 Total Lock Swivel Casters

Bar Activated Hi-Lo Control allows adjustment of the table height while keeping both hands on the patient and both feet on the floor. Control is designed to allow easy access from anywhere around the table without disrupting the treatment.

The AM-BAX3000 and AM-BAX5000 models are equipped with a motorized center section that allows positioning with the patient's full weight on the table. The center section motor controls are located on each side of the center section to allow adjustment without losing contact with the patient.