

AARON® SURGI-CENTER|PRO ELECTROSURGICAL GENERATOR

The Aaron Surgi-Center|PRO Electrosurgical Generator is intended for cutting, coagulation, ablation of tissue in general, gynecologic, orthopedic, ENT and urological procedures performed in an operating suite and procedure room. The unit was designed to provide an electrosurgery solution capable of most operating room demands in an affordable 200 watt generator. The unit offers 12 options of energy delivery for a variety of surgical application; 2 cut, 4 blend, 3 coagulation and 3 bipolar modalities.



Style # A2350

▼ Standard Equipment

Aaron® Surgi-Center|PRO unit
Hospital-grade power cord (120 VAC, 240 VAC)*
user's guide (MC-55-232-001)*

Dimensions

Height: 16.5cm (6.5 in.)
Width: 37.5cm (14.75 in.)
Depth: 46cm (18.1 in.)
Weight: <9.07 kg (<20 lbs)

Output Power

Cut I & II: 200 watts
Blend (1-4): 200 watts
Pinpoint Coagulation: 120 watts
Spray Coagulation: 120 watts
Gentle Coagulation: 120 watts
Macro Bipolar: 80 watts
Micro Bipolar: 80 watts
Standard Bipolar: 50 watts

Output Frequency

490 kHz ± 5

Line Voltage

100 – 240V ~ ± 10%

Line Frequency

50 – 60 Hz

Warranty: FOUR YEARS



* Not Pictured
† Sold Separately

★ RECOMMENDED ACCESSORIES †

To avoid incompatibility and unsafe operation, we recommend using the following Bovie® accessories with the Aaron® Surgi-Center|PRO:



BV-1253B



BV-1254B



A1255A



A827V



BV-IDS-CS



ESRE



A821



ESP1, ESP6, or ESPR2