

BrewerAccess High-Low EXAM TABLE



Abundant Storage Capacity

- Patented pass-through side drawers allow left or right access
- Eliminates right or left hand setup problems
- Additional spacious front drawer
- Convenient covered paper roll storage holds up to six rolls



Explore the Access High-Low family here. →



Chair



Table

Chair-To-Table Design

- Takes patient from consultation to examination
- 18" low height for easy and safe patient transfers
- Return to chair feature easily takes table back to chair position in one simple step (on power back model only)



Standard Patient Transfer Grab Bars

- Provides patients a stable and secure grab point to use when getting on and off table
- Conveniently swings out for lateral wheelchair transfers
- Easily removable



Pass-Through Work Surface

- Pulls out on either side
- Saves space by eliminating the need for additional stands or trays
- Perfect for charts, instruments, procedure preparation, or portable computer



Intuitive Foot Control

- Moisture resistant; Low voltage
- On / off switch located on the back of table allows table's foot control to be disabled and prevents patients and / or children from moving the table



Comfortable, Ergonomic Upholstery

- Seamless, medicinal stain resistant
- Available in 15 standard colors and UltraLeather™ plush upholstery
- Cal-133 fire code compliant upholstery available



Ergonomic Stirrups

- Standard, extendable, adjustable
- Conveniently stored in the table
- Manually extended through infinite positions and four lateral positions

Brewer

Furnishing solutions for safe practices.

BrewerAccess High-Low EXAM TABLE

Advanced Safety

- Safety grab bars for patient stability and positioning
- Return to chair feature on power back models
- Anti-tamper on/off switch prevents accidental movement
- Sensors prevent injuries during table movement
- 18" (45.7 cm) low height, 28" (71.1 cm) when flat
- Easy on / off for patients
- Backrest adjusts for easy patient positioning (pneumatic or power back available)

Patient Comfort

- Easy lateral wheelchair transfers
- Multi-density support upholstery
- Standard adjustable ergonomic stirrups
- Toe board / leg extension / stainless steel debris tray

Abundant Storage

- Pass-through work surface on side of table
- Patented pass-through side drawers accessible from right or left side of the table
- Front drawer: 19.8" x 19.2" x 3.8" (50.3 cm x 48.8 cm x 9.7 cm)
- Side drawers: 22.6" x 20.2" x 7" (57.4 cm x 51.3 cm x 17.8 cm)
- Covered storage for up to 6 additional paper rolls
- Paper roll storage capacity: 23" x 23" x 3.25" (58.4 cm x 58.4 cm x 8.3 cm)
- Slide through holder accepts 21" x 3" (53.3 cm x 7.6 cm) and 18" x 3" (45.7 cm x 7.6 cm) paper rolls

Certifications

All Brewer products conform to US and Canadian electrical and medical industry codes and regulations and meet or exceed industry requirements.

- UL 60601-1



Table Specifications

- 450 lb (204 kg) weight capacity
- Standard patient assist grab bars
- Standard adjustable stirrups
- Standard foot control
- Standard 3-year warranty
- Height: 18" - 37" (45.7 cm - 94 cm)
- Length: 57" extends to 78" (144.8 cm - 198.1 cm)
- Upholstery width: 28" (71.1 cm)
- Leg extension length: 21" (53.3 cm)
- Goes flat at 28" (71.1 cm)
- 115VAC, 60 Hz, 5.8 amps max
- Safe, low voltage actuators
- Safe, single power cord design, detachable
- Optional 5° pelvic tilt / 60-Watt front drawer heater
- Easy to clean, seamless upholstery
- 15 standard upholstery color options
- Cal-133 fire code compliant upholstery available in 4 colors
- Premium, plush upholstery available in 5 colors
- Custom colors also available

models

#6000	Access High-Low Exam Table with pneumatic back
#6001	Access High-Low Exam Table with pneumatic back, pelvic tilt & drawer heater
#6500	Access High-Low Exam Table with power back
#6501	Access High-Low Exam Table with power back, pelvic tilt, drawer heater & outlet

accessories

#98464	stainless steel treatment pan
#99502	one universal socket clamp
#99504*	deluxe articulating knee crutches
#100253	urology drain pan system
#100407	front drawer dividers (can order up to 7)
#100540	two support adapters for knee crutches (plant installed)
#100550	two support adapters for knee crutches (field installed)
#101800	electrical outlet option (6001 only)
#102031	Brewer Welch Allyn® IV bracket (plant installed)
#102032	Brewer Welch Allyn® IV bracket (field installed)

* Needs two #99502 and one of either #100540 or #100550

The Brewer Company is an ISO 13485:2003 certified company.

All Brewer Access™ High-Low Exam Tables come with a 3-year warranty.

The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

Call your local Brewer sales representative to learn about our complete product line.

The Brewer Company LLC
N88W13901 Main Street, Menomonee Falls, WI 53051
1.888.Brewer.1 | www.brewercompany.com

Brewer