



PL400H

Max Capacity 400Lb

- 400 Lb capacity hydraulic lift
- Six-Point spreader bar
- Manual Open/Close base legs
- 4" rear casters with locking brakes
- 3" front casters
- Tool free assembly
- Small & quick disassembly for transport
- Optional Upgrade to electric power
- Optional Digital weigh scale
- Low cost electronics standard
- Durable painted steel tube construction



Key dimensions of PL400H

NETWEIGHT		75Lb
A	base leg length	40"
B	base closed external width	26"
C	base open external width	37"
D	mast height	48"
E	base height	4.7"
F	min. lifting height	28"
G	max. lifting height	68"