

# PL500

Max Capacity 500Lb



- 500 Lb capacity electric power
- Increased lifting range and capacity
- Two-Point spreader bar standard
- Manual Open/Close base legs
- Manual emergency lowering
- 4" rear casters with locking brakes
- Dual 3" front casters w/ bumpers
- Control box w/ detachable battery pack
- Optional dual battery / wall charger
- Optional Digital weigh scale
- Optional six-point spreader bar
- Optional walking arms
- Available with lower cost electronics

Key dimensions of PL500

NET WEIGHT	110Lb
A base legs length	46.8"
B base closed external width	26"
C base open external width	41"
D mast height	50"
E base height	4.7"
F min. lifting height	22.4"
G max. lifting height	75.4"

