

SA600

Max Capacity 600Lb

- 600 Lb capacity electric power
- Increased capacity in same footprint
- Manual Emergency Lowering standard
- Manual base Open/Close
- 4" rear casters with locking brakes
- Dual 3" front casters with bumpers
- Control box w/ detachable battery pack
- Extra wide cushioned knee pad
- Extra wide removable foot rest
- Optional dual battery / wall charger
- Optional calf & buttock support strap
- Available with lower cost electronics
- Use with Standing or Toileting slings



Key dimensions of SA600

NET WEIGHT	139Lb
A base legs length	40"
B internal width	19"
C external width	26"
D expandable internal width	38"
E base height	4.7"
F min. lifting height	33"
G max. lifting height	68"

