

STA400

Max Capacity 400Lb



- 400 Lb Capacity
- Folding seats
- Tapered front legs for access
- Optional seat Locks
- 4" rear casters with locking brakes
- 3" front casters
- Durable, molded knee pads
- Rubber foot rest
- Grab bar for user assisted transfer
- Hooks for optional side support strap
- "S" model available with low seat

Key dimensions of STA400

NET WEIGHT		57Lb
A	min. base inner width	13.7"
B	max. base outer width	20"
C	fixed base outer back width	25"
D	seat width	20"
E	overall length	31"
F	seat height	26.3"
G	overall height	43"

