

VisionAire™ 2, 3 & 5

Compact Stationary Oxygen Concentrators

VisionAire™ 2 — For Pediatric Oxygen Needs

VisionAire 2 supports pediatric care with two types of precision flowmeters. The standard flowmeter has 1/8 LPM increments ranging from 1/8 LPM up to 2 LPM with CCM markings for easy conversion, and is suitable for adult use where lower power consumption and cooler operating temperatures are valued. The alternative flowmeter has 1/10 increments ranging from .10 LPM up to 1.0 LPM.

Lightest, Quietest, and Most Power-Efficient at 3 & 5 LPM

VisionAire 3 and 5 are the most power-efficient stationary oxygen concentrators on the market. Offering up to 3 LPM and 5 LPM, respectively, the VisionAire 3 and 5 operate at nearly silent decibel levels providing an added layer of patient comfort at home.

Patients and providers have long established their objectives concerning in-home oxygen. Only with the VisionAire™ series compact concentrators, including VisionAire 3 & 5, are all these objectives met.

Cost-Effective, Maintenance-Free Units

You can experience the added assurance of maintenance-free oxygen concentrators, fully backed by the expert in oxygen-generating technology. This helps ensure that you maximize the potential of your equipment inventory while greatly reducing or eliminating any unnecessary down time.

The technology of these filter-free units offers AirSep's time-proven expertise in a smaller cabinet that will please patients and their physicians, and meet your business requirements. Save time and money with the VisionAire Family!

With Safety in Mind

Important safety features include a backlined flowmeter for viewing accuracy and proper setting, a locking flow, and overdraw protection. Each concentrator is equipped with a resettable circuit breaker, a dual-function compliance/hour meter, and well-defined system alarms.

Features

- Most power-efficient 3 LPM (175 W) & 5 LPM (290 W) units on the market
- Compact, lightweight cabinet design — 30 lb (13.6 kg)
- Near silent operation with a low sound level of just 40 dBA
- Maintenance-free
- VisionAire 2 equipped with oxygen monitor. Optional for VisionAire 3 & 5
- Optional air outlet for nebulizer treatments for VisionAire 3 & 5



Call 1.800.482.2473 to order!
www.CAIREmedical.com • www.AirSep.com

VisionAire™ 2, 3 & 5

Compact Stationary Oxygen Concentrators



Product Specifications

| | VisionAire 2 | VisionAire 3 | VisionAire 5 |
|-----------------------|--|---|--|
| Oxygen Concentration* | Up to 1 LPM at 95.0% +0.5%/-3%; 2 LPM at 90% +5.5%/-3% | 0.5-3 liters per minute at 90% + 5.5/-3% | 1-5 liters per minute at 90% + 5.5/-3% |
| Dimensions In (cm) | 20.8 x 14.1 x 11.5 (52.8 x 35.8 x 29.2) | 20.8 x 14.1 x 11.5 (52.8 x 35.8 x 29.2) | 20.8 x 14.1 x 11.5 (52.8 x 35.8 x 29.2) |
| Weight | 30 lb; shipping weight - 37 lb (13.6 kg; shipping weight - 16.7 kg) | 30 lb; shipping weight - 37 lb (13.6 kg; shipping weight - 16.7 kg) | 30 lb; shipping weight - 37 lb (13.6 kg; shipping weight - 16.7 kg) |
| Electrical | 115 V, 60 Hz, 200 watts 230 V, 50 Hz, 200 watts 230 V, 60 Hz, 200 watts | 100 V, 50/60 Hz, 175 watts 115 V, 60 Hz, 175 watts 230 V, 50 Hz, 175 watts 230 V, 60 Hz, 175 watts | 115 V, 60 Hz, 3 amps, 290 watts 220-240 V, 50 Hz, 1.5 amps, 290 watts 230 V, 60 Hz, 1.5 amps, 290 watts |
| Sound Level | 40 dBA | 40 dBA | 40 dBA |
| Alarms | Power failure, High and low pressure, High temperature, No flow, Low concentration | Power failure, High and low pressure, High temperature, No flow, Low concentration (with optional oxygen monitor) | Power failure, High and low pressure High temperature, No flow, Low concentration (with optional oxygen monitor) |
| Relative Humidity | Up to 95% (non-condensing) | Up to 95% (non-condensing) | Up to 95% (non-condensing) |

*Based on an atmospheric pressure of 14.7 psia (101 kPa) at 70° F (21 °C). Operating unit outside of operational temperature range can affect performance

WARRANTY: 3-Year Standard Warranty. Other warranty options available.

INTERNATIONAL APPROVALS: CE 0459; ISO 9001:2008 and ISO 13485:2003-certified.

Ordering Information

| Description | Catalog Number** | Description | Catalog Number** |
|--|------------------|---|------------------|
| VisionAire 2 with oxygen monitor (115 V, 60 Hz unit) | AS096-8 | VisionAire 5 (115 V, 60 Hz unit) | AS098-1 |
| VisionAire 3 (230 V, 50 Hz unit) | AS096-1 | VisionAire 5 (220-240 V, 50 Hz unit) | AS098-2 |
| VisionAire 3 with oxygen monitor (230 V, 50 Hz unit) | AS096-2 | VisionAire 5 (230 V, 60 Hz unit) | AS098-3 |
| VisionAire 3 with oxygen monitor; UK-style power cord (230 V, 60 Hz unit) | AS096-3 | VisionAire 5 with oxygen monitor (115 V, 60 Hz unit) | AS098-4 |
| VisionAire 3 with oxygen monitor; Japan-style power cord (100 V, 50/60 Hz unit) | AS096-4 | VisionAire 5 with oxygen monitor (220-240 V, 50 Hz unit) | AS098-5 |
| VisionAire 3 with oxygen monitor (230 V, 60 Hz unit) | AS096-5 | VisionAire 5 with oxygen monitor (230 V, 60 Hz unit) | AS098-6 |
| VisionAire 3 (115 V, 60 Hz unit) | AS096-6 | VisionAire 5 with air outlet (115 V, 60 Hz unit) | AS098-707 |
| VisionAire 3 with oxygen monitor (115 V, 60 Hz unit) | AS096-7 | VisionAire 5 with air outlet (220-240 V, 50 Hz unit) | AS098-708 |
| | | VisionAire 5 with air outlet (230 V, 60 Hz unit) | AS098-709 |
| | | VisionAire 5 with oxygen monitor and air outlet (115 V, 60 Hz unit) | AS098-710 |
| | | VisionAire 5 with oxygen monitor and air outlet (220-240 V, 50 Hz unit) | AS098-711 |
| | | VisionAire 5 with oxygen monitor and air outlet (230 V, 60 Hz unit) | AS098-712 |

**Other model units available. Catalog Numbers reflect US English Part Numbers, for language-specific part numbers, please contact Customer Service.

CAIRE, Inc.
2200 Airport Industrial Drive, Suite 500
Ball Ground, GA 30107
Ph 770-721-7759 • Toll-Free 1-800-482-2473
Fax 888-932-2473

Chart BioMedical, Ltd.
Unit 6, Ashville Way
Wokingham RG41 2PL, United Kingdom
Ph +44(0) 1189 367060 • Fax +44(0) 1189 799245