

BLOOD DRAWING CHAIRS



Examine Real Quality™



1 Introduction

2—3 MD Series

4—5 Value Series

6—7 Standard Lab Series

8—10 Lab X Series

11 Bariatric Series

12 Specialty Items

13 Vinyl Color Selector

Product Index by Model Number

674	5
675	7
6000B-F	11
6000B-FD	11
6000X-F	10
6000X-FD	10
6010-F	9
6011-F	9
6011-V	5
6011-P	7
6020-F	8
6022-F	9
6022-V	5
6022-P	7
6029-V	4
6029-P	6
6060-V	4
6060-P	7
6077	12
6088	12
6099-F	10
6099-FD	10
64929-F	3
64929-BF	3
64950-BF	3
64950-F	2

EXAMINE REAL QUALITY

■ **QUALITY.** Even after over 25 years of manufacturing, Clinton never takes quality lightly. In fact, we invite examination. You'll find welded frames, reinforced stress points, safety-minded adjustments and tough, premium vinyl upholstery. That's real quality and that's what Clinton equipment is all about.

■ **INNOVATION.** Some would be satisfied just selling good products, but not Clinton. We take good products and make them better, stronger and easier to use with innovative ideas and state of the art precision manufacturing. It's not the easy way to make equipment, but it's the right way.

■ **CUSTOMER SERVICE.** You talk, we'll listen. It's really pretty simple. When you call Clinton, you talk to real people, not machines. Clinton customer service is second to none, and we intend to keep it that way.



MD SERIES

FEATURES

- Contemporary all-welded steel frame
- Adjustable height and depth flip arm
- Armrest locks in place when in use, then flips out of the way for easy exit
- Comfortable 2" padding throughout
- Easy-clean, antimicrobial, heavy weight, vinyl upholstery
- White powder-coated frame
- 2 knobs lock chrome armrest into position
- Non-marking glides
- Adjustable arm width
- 300 lbs. capacity

	Seat	Height
64950-F	17" x 16"	20"
	Width	Depth
Overall	32"	26"
	Arm Width Adjustment	Arm Height Adjustment
Arm	19"—24"	28"—33"

Blood Drawing Chair

- Padded and upholstered seat and armrests
- Adjustable arms
- White chair frame



Arm flips clear for easy seating and exit



Two knobs securely lock armrest into position for safety





	Seat	Height
64929-F	17" x 16"	20"
	Width	Depth
Overall	42"	26"
	Arm Width Adjustment	Arm Height Adjustment
Arm	19"—24"	28"—33"

Blood Drawing Chair with Drawer

- Padded and upholstered seat and armrests
- Adjustable arms
- White chair frame
- Easy-clean, one-piece plastic drawer



	Seat	Height
64929-BF	17" x 16"	20"
	Width	Depth
Overall	42"	26"
	Arm Width Adjustment	Arm Height Adjustment
Arm	19"—24"	28"—33"

Blood Drawing Chair with Drawer

- Padded and upholstered seat, back and armrests
- Adjustable arms
- White chair frame
- Easy-clean, one-piece plastic drawer



	Seat	Height
64950-BF	17" x 16"	20"
	Width	Depth
Overall	32"	26"
	Arm Width Adjustment	Arm Height Adjustment
Arm	19"—24"	28"—33"

Blood Drawing Chair

- Padded and upholstered seat, back and armrests
- Adjustable arms
- White chair frame

VALUE SERIES

FEATURES

- 1" square, all-welded steel frame
- 360° non-locking swing arm
- Height adjustable armrests
- All laminate arms with tee-molded edge bumpers
- Easy-clean, one-piece plastic seat
- Optional padded armrests
- Neutral Grey powder-coated frame
- 2 knobs lock chrome armrest into position
- Adjustable levelers
- 300 lbs. capacity



	Seat	Height
6029-V	17" x 16"	20"
	Width	Depth
Overall	43"	26"
	Height Adjustment	
Arms	27"—33"	

Blood Drawing Chair with Drawer

- One-piece plastic seat
- Adjustable arms
- Laminate pivot armrest
- Easy-clean, one-piece drawer

	Seat	Height
6060-V	17" x 16"	20"
	Width	Depth
Overall	34"	26"
	Height Adjustment	
Arms	27"—33"	

Blood Drawing Chair

- One-piece plastic seat
- Adjustable arms
- Laminate pivot armrest





	Seat	Height
6011-V	17" x 16"	26"
	Width	Depth
Overall	34"	26"
	Height Adjustment	
Arms	34"—40"	

Extra-Tall Blood Drawing Chair

- Reduces bending for staff
- One-piece plastic seat
- Adjustable arms
- Laminate pivot armrest
- Patient foot rest
- Makes it easy for patient to sit down and exit



	Seat	Height
6022-V	17" x 16"	26"
	Width	Depth
Overall	43"	26"
	Height Adjustment	
Arms	34"—40"	

Extra-Tall Blood Drawing Chair with Drawer

- Reduces bending for staff
- One-piece plastic seat
- Adjustable arms
- Laminate pivot armrest
- Patient foot rest
- Easy-clean, one-piece plastic drawer
- Makes it easy for patient to sit down and exit



674	Option
------------	--------

Optional Padded Armrest

- Same as standard laminate armrest with the added comfort of upholstery
- Please specify color when ordering

STANDARD LAB SERIES

FEATURES

- 1 1/4" square, heavy duty, all-welded, tubular steel frame
- Adjustable height and depth flip arm
- Armrest locks in place when in use, then flips out of the way for easy exit
- Neutral Grey powder-coated frame
- 2 knobs lock chrome armrest into position
- Adjustable foot levelers
- Complete series of chairs to choose from
- 300 lbs. capacity



One-piece flip-arm with edge-bumper stripping



2 knobs lock chrome armrest into position

	Seat	Height
6029-P	17" x 16"	20"
	Width	Depth
Overall	43"	26"
	Height Adjustment	
Arms	27"—33"	

Blood Drawing Chair with Flip Arm and Drawer

- Easy-clean, one-piece plastic seat
- Adjustable laminate flip-arm
- Easy-clean one piece drawer
- Optional padded arms





	Seat	Height
6060-P	17" x 16"	20"
	Width	Depth
Overall	34"	26"
	Height Adjustment	
Arms	27"—33"	

Blood Drawing Chair with Flip Arm

- Easy-clean, one-piece plastic seat
- Adjustable laminate flip-arm
- Optional padded arms

Arm in up position



	Seat	Height
6011-P	17" x 16"	26"
	Width	Depth
Overall	34"	26"
	Height Adjustment	
Arms	34"—40"	

Extra-Tall Blood Drawing Chair with Flip-Arm

- Easy-clean, one-piece plastic seat
- Adjustable laminate flip-arm
- Optional padded arms



	Seat	Height
6022-P	17" x 16"	26"
	Width	Depth
Overall	43"	26"
	Height Adjustment	
Arms	34"—40"	

Extra-Tall Blood Drawing Chair with Flip-Arm and Drawer

- Easy-clean, one-piece plastic seat
- Adjustable laminate flip-arm
- Easy-clean, one-piece drawer
- Optional padded arms



675 Option

Optional Padded Flip-Arm for Standard Lab Series

- 2" padded armrests
- Please specify color when ordering

LAB X SERIES

FEATURES

- 1 1/4" square, heavy duty, all-welded, tubular steel frame
- Adjustable height and depth flip arm
- Armrest locks in place when in use, then flips out of the way for easy exit
- Comfortable 2" padding throughout
- Neutral Grey powder-coated frame
- 2 knobs lock chrome armrest into position
- Adjustable foot levelers
- Easy-clean durable antimicrobial vinyl
- Complete series of chairs to choose from
- 400 lbs. capacity



Padded, adjustable flip-arm



2 knobs lock chrome armrest into position

	Seat	Height
6020-F	17" x 16"	20"
	Width	Depth
Overall	43"	26"
	Height Adjustment	
Arms	27"—33"	

Blood Drawing Chair with Flip Arm and Drawer

- Adjustable, padded and upholstered flip-arm
- Easy-clean, one-piece drawer





	Seat	Height
6010-F	17" x 16"	20"
	Width	Depth
Overall	34"	26"
	Height Adjustment	
Arms	27"—33"	

Blood Drawing Chair with Flip Arm

- Adjustable, padded and upholstered flip-arm



	Seat	Height
6022-F	17" x 16"	26"
	Width	Depth
Overall	43"	26"
	Height Adjustment	
Arms	34"—40"	

Extra-Tall Blood Drawing Chair with Flip Arm and Drawer

- Adjustable, padded and upholstered flip-arm
- Easy-clean, one-piece drawer
- Extra height reduces bending for staff
- Easy for patient to sit and exit chair
- Patient footrest

	Seat	Height
6011-F	17" x 16"	26"
	Width	Depth
Overall	34"	26"
	Height Adjustment	
Arms	34"—40"	

Extra-Tall Blood Drawing Chair with Flip Arm

- Adjustable, padded and upholstered flip-arm
- Extra height reduces bending for staff
- Easy for patient to sit and exit chair
- Patient footrest





	Seat	Height
6000X-F	35" x 18"	20"
	Width	Depth
Overall	48"	28"
	Height Adjustment	
Arms	29"—34"	



One arm in up position

Extra-Wide Blood Drawing Chair with Flip-Arms

- Padded and upholstered seat, back and armrests
- 2 adjustable flip-arms
- Accommodates large patient or child with parent up to 400 lbs. combined weight



	Seat	Height
6099-F	35" x 18"	26"
	Width	Depth
Overall	48"	28"
	Height Adjustment	
Arms	34"—40"	



One arm in up position

Extra-Tall & Wide Blood Drawing Chair with Flip-Arms

- Padded and upholstered seat, back and armrests
- 2 adjustable flip-arms
- Extra height reduces bending for staff
- Accommodates large patient or child with parent up to 400 lbs. combined weight



	Seat	Height
6000X-FD	35" x 18"	20"
	Width	Depth
Overall	58"	28"
	Height Adjustment	
Arms	29"—34"	

Extra-Wide Blood Drawing Chair with Drawer and Flip-Arm

- Padded and upholstered seat, back and armrests
- 1 adjustable flip-arm
- Easy-clean, one-piece drawer
- Accommodates large patient or child with parent up to 400 lbs. combined weight



	Seat	Height
6099-FD	35" x 18"	26"
	Width	Depth
Overall	58"	28"
	Height Adjustment	
Arms	34"—40"	

Extra-Tall & Wide Blood Drawing Chair with Drawer and Flip-Arm

- Padded and upholstered seat, back and armrests
- 1 adjustable flip-arm
- Easy-clean, one-piece drawer
- Extra height reduces bending for staff
- Accommodates large patient or child with parent up to 400 lbs. combined weight

BARIATRIC SERIES

FEATURES

- Seat has additional depth and an extra center support for added strength
- Extra dimensions make it easier for patient and technician
- Adjustable height and depth flip arms
- Armrests lock in place when in use, then flips out of the way for easy exit
- Comfortable 2" padding throughout
- Easy-clean, antimicrobial, heavy weight, vinyl upholstery
- Neutral Grey powder-coated frame
- 2 knobs lock chrome armrest into position
- Adjustable foot levelers
- Chairs may be special ordered with straight armrests instead of flip arms (677 Option)
- Accommodates patients up to 700 lbs.



	Seat	Height
6000B-F	35" x 21"	20"
	Width	Depth
Overall	48"	32"
	Height Adjustment	
Arms	29"—34"	

Bariatric Blood Drawing Chair

- Accommodates larger patient or child with parent up to 700 lbs. combined weight
- Padded and upholstered seat, back and armrests
- 2 adjustable flip-arms

	Seat	Height
6000B-FD	35" x 21"	20"
	Width	Depth
Overall	58"	32"
	Height Adjustment	
Arms	29"—34"	

Bariatric Blood Drawing Chair with Drawer and Flip-Arm

- Accommodates larger patient or child with parent up to 700 lbs. combined weight
- Padded and upholstered seat, back and armrests
- 1 adjustable flip-arm
- Easy-clean, one-piece drawer



SPECIALTY ITEMS



Conveniently stores against wall out of the way of foot traffic.

	Depth Adjustment	Height Adjustment
6088	24"—36"	26"—39"
	Width	Length
Armrests	6"	24"

Wheelchair Blood Drawing Station

- Wall-mounted for secure dependable performance
- Adjustable height and depth
- Armrest locks in place when in use, then flips out of the way for easy exit
- Comfortable 2" padding
- Easy-clean, antimicrobial, heavy-weight, vinyl upholstery
- 2 knobs lock chrome armrest into position
- Plastic work tray (optional)

	Width	Length	Height
6077	20 1/2"	39"	39 1/2"
	Width	Length	
Top	20"	32"	

Infant Blood Drawing Station

- All-welded stainless steel frame
- Top has 2 rims
- 2 safety straps
- Plastic work tray
- Removable middle padded section



VINYL COLOR SELECTOR

ABOUT OUR UPHOLSTERY

Clinton uses only premium, woven knit-backed upholstery that resists stretching and bunching associated with ordinary backed vinyl upholstery. Plus, Clinton upholstery is topcoated with Perma-Block-3 to resist stains, scuffs and abrasions.



Slate Blue

3SB



Desert Tan

3DT



Tomato

3TM



China Green

3CG



Yellow

3YW



Black

3BK



Gunmetal

3GM



Dove Gray

3DG



Cream

3CM



Clamshell

3CS



Wedgewood

3WW



Royal Blue

3RB



Emerald

3EM




Mulberry

3MB



Purplegray

3PG



Burgundy

3BG

PLEASE NOTE:

Due to the limitations of the printing process, actual colors may vary.

Clinton Industries reserves the right to improve or change our products without prior notification.

See the entire Clinton product line. Call Clinton today for any of our product catalogs.



Medical Tables & Cabinets



Stools & Medical Accessories



Pediatric Medical Equipment



Blood Drawing Chairs



Physical Therapy Equipment



Clinton Industries 3-Year Limited Warranty

Products backed by this warranty are warranted against defects in workmanship for three years from date of purchase. This warranty does not include electrical products. Contact Clinton industries for full warranty details.



Examine Real Quality™

1140 East Edison St., York, PA 17402

800-441-9131 ■ 717-848-3519

Toll-Free Fax 866-522-1233

Fax 717-843-5871

www.clinton-ind.com