



**Wheelchair Blood Drawing Station**

- Wall-mounted for secure, dependable performance
- Adjustable height and depth
- Armrest positions in place when in use, then flips out of the way for easy exit
- Comfortable 2" padding
- Easy-clean, heavyweight, vinyl upholstery
- Two knobs position chrome armrest

|             |                  |                   |
|-------------|------------------|-------------------|
| <b>6088</b> | depth adjustment | height adjustment |
| overall     | 24" – 36"        | 26" – 39"         |
|             | 61 cm – 91.44 cm | 66 cm – 99 cm     |
| armrests    | width            | length            |
|             | 6"               | 24"               |
|             | 15.24 cm         | 61 cm             |



Conveniently stores against wall out of the way of foot traffic.

**Half Round, Stationary, Padded Phlebotomy Stand**

- Vinyl upholstered, padded surface
- Easy-clean, 10<sup>3</sup>/<sub>4</sub>" L x 8<sup>1</sup>/<sub>4</sub>" W x 3<sup>1</sup>/<sub>2</sub>" D plastic bin (28 cm L x 21 cm W x 9 cm D)
- Chrome plated steel pole with knob-lock height adjustment
- 5-leg, cast aluminum base with black, non-marking, nylon glides

|               |                    |              |
|---------------|--------------------|--------------|
| <b>6940</b>   | depth              | radius width |
| top           | 22"                | 11"          |
|               | 55.88 cm           | 27.94 cm     |
| base diameter | 23"                |              |
|               | 58.42 cm           |              |
| height range  | 30" – 38"          |              |
|               | 76.2 cm – 96.52 cm |              |



(Supplies shown on all models not included)

**Half Round, Stationary, ClintonClean™ Phlebotomy Stand**

- ClintonClean™ solid plastic top
- Easy-clean, 11" L x 8<sup>1</sup>/<sub>4</sub>" W x 3<sup>1</sup>/<sub>2</sub>" D plastic bin (28 cm L x 21 cm W x 9 cm D)
- Chrome plated steel pole with knob-lock adjustment
- 5-leg, cast aluminum base with black, non-marking, nylon glides

|               |                    |                                  |
|---------------|--------------------|----------------------------------|
| <b>6950</b>   | depth              | radius width                     |
| top           | 22"                | 11 <sup>1</sup> / <sub>2</sub> " |
|               | 55.88 cm           | 29.21 cm                         |
| base diameter | 23"                |                                  |
|               | 58.42 cm           |                                  |
| height range  | 30" – 38"          |                                  |
|               | 76.2 cm – 96.52 cm |                                  |



6940 Demo

6950 Demo

