

FLAT TOP TREATMENT TABLES





Classic Design. Premium Quality. Exceptionally Priced.

Not all treatment tables are created equal — or priced equally.

Clinton Flat Top Tables use materials like solid hardwood legs, premium 1 1/8" panel legs or a sleek, powder-coated steel frame. These quality materials aren't inexpensive, but Clinton's classic designs demand nothing less. As a matter of fact, some Clinton treatment tables have been in service for more than 25 years. We credit much of our durability to our use of steel fastener

construction instead of wooden dowels, which can weaken over time. All of our tables feature a comfortable 2" firm foam padded top in your choice of 20 colors of premium, woven-back upholstery. All the quality, all the selection — and all at some of the best prices in the business, PLUS a 5-year warranty! Examine our flat top tables, and you'll find why Clinton flat out beats the competition in quality, selection and price.

CLASSIC SERIES  
TREATMENT TABLES



Classic Series, Wood, Treatment Tables feature:

- Solid hardwood legs
- H-brace construction
- Legs have pre-inserted carriage bolts for fast easy assembly
- Seamless, rounded corner top
- 2" of firm foam padding (5 cm)
- Available in Natural and Dark Cherry finish
- Choice of 20 upholstery colors
- 5-year, limited warranty
- Fast, easy assembly with basic tools
- 350 lbs. load capacity under normal use (158.75 kg)

Straight Line Treatment Table

100-24	length	width	height
	72"	24"	31"
	182.88 cm	60.95 cm	78.73 cm

100-27	27" width model (68.58 cm)
--------	----------------------------

100-30	30" width model (76.19 cm)
--------	----------------------------

Straight Line Treatment Table with Full Shelf

200-24	length	width	height
	72"	24"	31"
	182.88 cm	60.95 cm	78.73 cm

200-27	27" width model (68.58 cm)
--------	----------------------------

200-30	30" width model (76.19 cm)
--------	----------------------------

Straight Line Treatment Table with Shelving

300-24	length	width	height
	72"	24"	31"
	182.88 cm	60.95 cm	78.73 cm

300-27	27" width model (68.58 cm)
--------	----------------------------

300-30	30" width model (76.19 cm)
--------	----------------------------

Straight Line Treatment Table with Shelf  
& Two Drawers

1300-27	length	width	height
	72"	27"	31"
	182.88 cm	68.58 cm	78.73 cm

1300-30	30" width model (76.19 cm)
---------	----------------------------

Alpha S-Series, Steel, Treatment Tables feature:

- Gray, powder-coated, all-welded, steel frame
- Adjustable height levelers
- Pre-threaded hardware for easy assembly
- Paper dispenser included
- 2" of firm foam padding (5 cm)
- Choice of 20 upholstery colors
- Traditional square corner top
- 5-year, limited warranty
- Fast, easy assembly with basic tools
- 500 lbs. load capacity under normal use (226.79 kg)

Straight Line Treatment Table

3100-27	length	width	height
	72"	27"	31"
	182.88 cm	68.58 cm	78.73 cm

3100-30	30" width model (76.19 cm)
---------	----------------------------

Straight Line Treatment Table with Full Shelf

3200-27	length	width	height
	72"	27"	31"
	182.88 cm	68.58 cm	78.73 cm

3200-30	30" width model (76.19 cm)
---------	----------------------------

Straight Line Treatment Table with Shelving

3300-27	length	width	height
	72"	27"	31"
	182.88 cm	68.58 cm	78.73 cm

3300-30	30" width model (76.19 cm)
---------	----------------------------

Straight Line Treatment Table with Shelf  
& Two Drawers

3130-27	length	width	height
	72"	27"	31"
	182.88 cm	68.58 cm	78.73 cm

3130-30	30" width model (76.19 cm)
---------	----------------------------



ALPHA S-SERIES  
TREATMENT TABLES



## STYLE LINE TREATMENT TABLES

### *Style Line, Panel Leg, Treatment Tables* feature:

- Steel reinforced, pre-assembled plywood top frame
- 1 1/8" panel legs with nylon glides (2.85 cm)
- Vertical spanner for extra support
- Bolt-through, pre-threaded hardware
- Easy-clean laminate surfaces
- Chip resistant PVC edges
- Seamless, rounded corner top
- 2" of firm foam padding (5 cm)
- Available in your choice of select laminate colors
- Choice of 20 upholstery colors
- 5-year, limited warranty
- Fast, easy assembly with basic tools
- 350 lbs. load capacity under normal use (158.75 kg)



9300



9313

### Straight Line Treatment Table

9100-27	length	width	height
	72"	27"	31"
	182.88 cm	68.58 cm	78.73 cm

9100-30 30" width model (76.19 cm)

### Straight Line Treatment Table with Full Shelf (pictured on cover)

9200-27	length	width	height
	72"	27"	31"
	182.88 cm	68.58 cm	78.73 cm

9200-30 30" width model (76.19 cm)

### Straight Line Treatment Table with Shelving

9300-27	length	width	height
	72"	27"	31"
	182.88 cm	68.58 cm	78.73 cm

9300-30 30" width model (76.19 cm)

### Straight Line Treatment Table with Shelf and Two Drawers

9313-27	length	width	height
	72"	27"	31"
	182.88 cm	68.58 cm	78.73 cm

9313-30 30" width model (76.19 cm)



9100





## UPHOLSTERY

Clinton uses tough, woven, knit-backed, vinyl upholstery that resists bunching and stretching. Its durable surface can be sanitized quickly and easily with your specific disinfectant and a soft cloth.

*(See detailed cleaning instructions with product.)*

In an effort to be current on changing color trends and to better provide our customers with the most up-to-date palette, please visit our web site at [www.clinton-ind.com](http://www.clinton-ind.com).

A selection of Clinton upholstery colors available for each specific product can be found on every product page by clicking the "Color Chart" tab.

## LAMINATE

In an effort to be current on changing color trends and to better provide our customers with the most up-to-date palette, please visit our web site at [www.clinton-ind.com](http://www.clinton-ind.com).

A selection of Clinton laminate available for each specific product can be found on every product page by clicking the "Color Chart" tab.

## DOOR & DRAWER PULLS

**Door and drawer pulls** can coordinate with solid color base laminates. Woodtone base colors look best with silver or black wire pulls. If no pulls are specified, Gray pulls will be used on solid color fronts and Silver Wire pulls will be used on woodtones.

In an effort to be current on changing trends and to better provide our customers with the most up-to-date hardware, please visit our web site at [www.clinton-ind.com](http://www.clinton-ind.com).

A selection of Clinton door and drawer pulls available for each specific product can be found on every product page by clicking the "Color Chart" tab.

[Specs.](#)[Video](#)[Options](#)[Color Chart](#)[Catalog](#)[Related Products](#)

800-441-9131 • 717-848-3519

Toll-Free Fax 866-522-1233 • Fax 717-843-5871

[www.clinton-ind.com](http://www.clinton-ind.com)

