




Bend-A-Beam® ABC® Handpieces


- Reduced blood loss
- Effective on high impedance tissue such as bone
- Reduced smoke plum for better visualization



Bend-A-Beam® Malleable ABC® Handpieces

High-quality disposable Argon Beam Coagulation handpieces are designed to deliver Argon Beam Coagulation to remote surgical sites.

- Bendable up to 60°
- Available in 3 inch, 6 inch and 9 inch versions, insulated and malleable.



For use with the System 7550™
Electrosurgical Generator with Argon
Beam Coagulation

Not intended for laparoscopic use.
Refer to the Instructions for Use for specifics.



Ordering Information

Catalog Number	Product Description	Packaging
134003	3" (7.6cm) Bend-a-Beam® Single Function Malleable ABC® Hand Control Handpiece with 10' (3.05m) Cord. Allows Use of Monopolar ABC® Coagulation Only. Single Use, Sterile	10/case 1/pkg
134006	6" (15.2cm) Bend-a-Beam® Single Function Malleable ABC® Hand Control Handpiece with 10' (3.05m) Cord. Allows Use of Monopolar ABC® Coagulation Only. Single Use, Sterile	10/case 1/pkg
134009	9" (22.8cm) Bend-a-Beam® Single Function Malleable ABC® Hand Control Handpiece with 10' (3.05m) Cord. Allows Use of Monopolar ABC® Coagulation Only. Single Use, Sterile	10/case 1/pkg

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



525 French Road
Utica, NY 13502
(315) 797-8375
Customer Service: 1-800-448-6506
FAX: 1-800-438-3051
International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235
email: info@mail.conmed.com
www.conmed.com

Please contact your local ConMed Territory Manager for information on our full line of Argon Beam Coagulation delivery devices and accredited ABC® Technology Continuing Education Programs.