Model 2441

Mechanical Pediatric and Infant Scale



Perfect for hospitals, pediatrician's offices, clinics, or the home nursery, Detecto's model 2441 pediatric and infant scale features a 16 kg capacity with 100 g increments. These sturdy mechanical scales are engineered with the active baby in mind and are built with the uncompromising quality you expect from Detecto.

The white metal tray holds infants comfortably, while the beam is easy to read with blue graduations and numerals on a white background. This scale also comes in 20 lb x 2 oz and 10 kg x 50 g models.

MODEL 2441 FEATURES

- Springless, balance beam
- Easy-to-read, blue graduations on white background
- Standard weight also available

Economical, compact beam scale

MODEL 2441 SPECIFICATIONS	
Capacity	16 kg x 100 g
Shipping Weight	21 lb / 10 kg
Tray Size	11" W x 19" L x 2 1/2" D 280 mm W x 483 mm L x 64 mm D
Overall Construction	10 1/2" H x 21" W 19" D 267 mm H x 533 mm W x 483 mm D
Metal Tray Construction	White Epoxy

A Division of Cardinal Scale Manufacturing Co.

203 E. Daugberty, Webb City, MO 64870 USA
Pb: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001

www.detectoscale.com

BUILDING A BETTER WEIGH