





D

E

E

Detecto

Model

6337

EYE-LEVEL DIGITAL HEALTH CARE SCALES

DETECTO MODELS 6337/6339

Detecto's ProMed® scales are the right choice for home, health, and clinical office weighing. These all-USA-manufactured scales are designed for years of dependable and accurate weighing, plus they are versatile enough to meet all clinical health requirements. With a capacity of up to 500 lbs, the 6337 and 6339 are designed to monitor weight in the home or medical office. Prescribed by dieticians and physicians, the inclusive Body Mass Index feature is easy to use for the monitoring of personal health/medical needs.

The economical CONVENIENCE of a USA-MADE **EYE-LEVEL DISPLAY** with **EXTRA-LARGE DIGITS**

The steel construction of Detecto's ProMed® scales ensure years of reliable weighing. The baked-on epoxy paint finish combines with a lift-off platform cover for easy cleaning to maintain a sanitary environment in the home or office. With the push of a button, readings appear in pounds or kilograms.

HEAVY SCALE BASE

Spacious weighing platform with heavy steel base for patient assurance during measurement.



MADE IN USA

Manufactured in the heart of America. Built with allsteel construction for longlasting weighing durability.



BODY MASS INDEX

Simply touch the BMI key and utilize the toggle arrows to adjust the patient height for BMI calculation.



VERSATILE POWER

Use either 6 "AA" batteries (not included) for efficient mobility or the AC adapter (optional).



HEIGHT ROD

Detecto's ProMed® model 6339 also offers an integral height rod.



Detecto

Model

6339

with

Height

Rod

PROMED



base as the standard in the medical weighing industry.

Sold By:

Heavy-duty 500lb capacity with medical grade accuracy to 0.2 lb.

Durable - long lasting - great investment.

6337 / 6339 SPECIFICATIONS **CAPACITY** 500 lb x 0.2 lb / 225 kg x 0.1 kg **HEIGHT ROD** 6337 None 6339 30" - 78" / 76 cm - 200 cm **POWER** 6 "AA" size Alkaline, Ni-Cad or NiMH batteries (not included) OR an optional 100 to 240 VAC 50/60Hz 12 VDC 1A wall plug-in UL/CSA listed AC power supply (not included) 10.75" W x 14.75" D / 27 cm W x 37.5 cm D **PLATFORM DISPLAY** Six-digit, seven-segment, 0.75" / 19 mm high LCD **KEYBOARD** Membrane type with 5 keys and 4 directional arrows CONSTRUCTION Painted steel **SCALE HEIGHT** 59" / 1.5 m SHIPPING WEIGHT 6337 34 lb / 15 kg 6339 37 lb / 16 kg

Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.







203 E. Daugherty, Webb City, MO 64870 USA Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001 www.Detecto.com

© Copyright 2009 Cardinal Scale Mfg. Co. • Printed in USA • CAR/00/0709/C228

PROMED

Detecto

Model

6337