

DIGITAL PEDIATRIC SCALES



**HEAVY-DUTY
SOLID, STEEL
SCALE BASE**
*for patient
assurance while
weighing.*

Easy-to-clean polystyrene trays cradle babies securely on these dependable pediatric scales.



***Digital speed
and accuracy***



**DETECTO
Models**

8435 | 8435KG

2 models to choose from

DETECTO DIGITAL PEDIATRIC SCALES




The Detecto brand name, known around the world, assures you of a reliable, accurate scale for pediatric weighing. The Detecto models 8435 and 8435KG offer you a choice of pounds and kilograms or kilograms. Lightweight and attractive, these models feature a heavy-duty, rust-resistant understructure and digital speed and accuracy. The polystyrene trays are a snap to clean and the scales feature compact footprints.

Detecto's model 8435 and 8435KG feature easy-to-use keypads and extra-large displays perfect for patient weighing and viewing by the attending nurse, physician, or home healthcare professional.



FEATURES:

- Pushbutton tare.
- Made in USA for convenient parts and service requirements.
- AC or battery power.
- Easy-to-clean plastic trays for baby comfort.
- Easy-to-read high-contrast LCD.
- Automatic weight hold feature for active babies, locks in weight.
- Heavy steel base for patient stability.
- Durable, long-lasting metal construction.



SPECIFICATIONS

CAPACITIES:	
8435	130 lb x 0.1 lb / 59 kg x 0.05 kg
8435KG	59 kg x 0.05 kg
OVERALL DIMENSIONS:	37.5" H x 22" W x 21" D / 95 cm H x 56 cm W x 53 cm D
TRAY DIMENSIONS:	4.5" H x 22" W x 14.75" D / 11 cm H x 56 cm W x 37 cm D
SHIPPING WEIGHT:	33 lb / 15 kg
CONSTRUCTION:	Baked-on, powder-painted steel
DISPLAY:	6-digit, 0.7"/18mm high-contrast LCD
POWER:	6 "AA" size alkaline batteries (not included) or optional AC adapter

Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.

 **DETECTO**[®]
A Division of Cardinal Scale Manufacturing Co.
203 E. Daugherty, Webb City, MO 64870 USA
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001
www.Detecto.com

Sold By:

ISO
REGISTERED 9001

