


HEALTH CARE
**STAINLESS
STEEL**
SCALES



Model 339S

Available in the world's standard,
the Eye-Level Mechanical Beam
Physician Scale and also in the
state-of-the-art, EMR-ready
Waist-High Digital Scale.

**YOUR CHOICE -
A WINNER EITHER WAY!**

**UNIQUELY DETECTO -
UNIQUELY STAINLESS STEEL**



Model 8437S



- **Hygiene** The time-proven benefits of *stainless steel* - cleans easily and completely.
- **Beauty** Shiny, attractive *stainless steel* - fits well in all modern decors.
- **Protection** Rust-resistant *stainless steel* protects the scale - especially useful in high humidity areas.

Stainless Steel MECHANICAL Health Care Scales

Enjoy the classic elegance of stainless steel in your Detecto mechanical weighbeam scale. A platform cover that lifts off easily combined with the stainless steel casing makes cleanup a breeze. Water damage from mopping and high humidity areas around swimming pools is a thing of the past with Detecto's stainless steel health care scales. Scale height is 59" (1.5 m) and the die-cast beam can be easily read from either side.



Model	Capacity		Height Rod 30 - 78 in 76 - 200 cm	Shipping Weight
	Avoirdupois	Metric		
439S	400 lb x 4 oz		•	38 lb / 17.2 kg
2391S		180 kg x 100 g	•	37 lb / 16.7 kg
437S	400 lb x 4 oz			35 lb / 15.8 kg
2371S		180 kg x 100 g		35 lb / 15.8 kg
339S	400 lb x 4 oz	175 kg x 100 g	•	38 lb / 17.2 kg
337S	400 lb x 4 oz	175 kg x 100 g		35 lb / 15.8 kg

Now, for the first time, you can combine the classic Detecto physician scale with the time-honored beauty and protection of stainless steel.

Stainless Steel DIGITAL Health Care Scales

With an ultra-modern design, Detecto's stainless steel digital health care scales offer an extra-tough finish which stands up to water damage from mopping around the base and high-humidity areas such as locker rooms. Enjoy the power of digital weighing with the convenience of "at-a-glance" weight readouts and an AC or battery-powered LCD display. With the push of a button, readings appear in pounds OR kilograms.

Model	Capacity 500 lb x 0.2 lb / 225 kg x 0.1 kg	Height Rod 30 - 78 in 76 - 200 cm	Shipping Weight
8430S	•	•	34 lb / 15.4 kg
8437S	•		31 lb / 14 kg



Detecto Scale reserves the right to improve, enhance or modify features and specifications without prior notice. All registered trademarks are the property of their respective owners.



SOLD BY:



Made in USA