



Mariner® WPS10DT

Wet Diaper Scale



**MEETING
NEONATAL
FACILITY
REQUIREMENTS,**

**RIGHT
DOWN TO THE
GRAM**

- 160 oz x 0.05 oz / 5,000 g x 1 g
- Tare Button to Remove Dry Diaper Weight and Measure Fluid Output
- 80-hour Rechargeable Battery Pack with Automatic Shut Off Feature
- 5/8-in Deep Stainless Steel Removable Tray Aids Sample Retention and Allows Easy Clean Up
- 1.2-in/30-mm High Backlit LCD Easily Readable in Dimly-Lit Nurseries
- Scale is Completely Submersible Underwater



DEPENDABLE ACCURACY for CLINICAL SETTINGS



DETECTO's
Mariner® WPS10DT wet
diaper scale meets all the hygiene
demands of any active medical environment.

EASY-TO-USE, PORTABLE, AND CONVENIENT

The backlit blue LCD display offers large 1-in/25-mm high bold digits, readily visible in even the dimmest nursery conditions. The raised buttons allow for tactile functionality and easy operation. The UNIT key toggles through six different weight capacity display options in both decimals and fractions. The tare function conveniently allows the weight of a dry diaper to be removed from the total weight of a soiled diaper. The convenient 80-hour rechargeable battery pack with automatic shut-off feature conserves energy and allows the scale to be moved where it is needed the most.



Removable stainless
steel tray makes
sanitation effortless.

HYGIENIC AND HEAVY-DUTY


The watertight AC plug and stainless steel construction allow for easy clean up, prevent corrosion, and reduce the risk of bacteria harboring in nooks and crannies. The 5/8-in deep removable stainless steel tray with raised edges holds items neatly on the weighing platform and makes clean up effortless. Both of these features are a necessity for messy nursery applications.

Backlit LCD is easily
visible, even in
darkened rooms
so babies won't
be disturbed.



Model	WPS10DT Wet Diaper Scales		Dimensions	Scale	
Capacity x Increments	10 lb x 0.1 oz	160 oz x 0.05 oz		Tray	7.5 in W x 9.8 in D x 2.3 in H / 190.5 mm W x 247.5 mm D x 57.0 mm H
	10 lb x 0.005 lb	160 oz x 1/8 oz			9.0 in W x 6.0 in D x 1.1 in H / 228.6 mm W x 152.4 mm D x 27.9 mm H
	10 lb x 1/8 oz	5,000 g x 1 g			
Zero	Established on power up routine and maintained by auto-zero circuitry		Tare		100% of full scale capacity
Display / Digits	1.2-in / 30-mm high backlit blue LCD, (5) 1-in / 25-mm high black LCD digits		Keypad Buttons		ON/OFF, UNITS/Backlight, SET, TARE
Feet	6 rubberized feet for traction		Enclosure Rating		IP67 submersible
Product Weight	4.9 lb / 2.2 kg		Operation Temperature		40° to 105°F (5° to 40°C)
Shipping Weight	7.1 lb / 3.2 kg		Storage Temperature		40° to 140°F (5° to 60°C)
Power	AC Adapter (included): 120 VAC input 9V, 100mA output - Battery (included): Ni-MH rechargeable 4.8 VDC 1800 mA		Please Note: The WPS10DT Mariner® is not commercial dishwasher safe.		
			Battery Life		80 hours of continuous use with a fully recharged battery pack
UPC Code	809161194203		Certification		NSF Approved

Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.

 **DETECTO**
A Division of Cardinal Scale Manufacturing Co.
203 E. Daugherty, Webb City, MO 64870 USA
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001
www.Detecto.com

Mobile: m.detecto.com

SOLD BY: