

Model DS8150

Digital Chair Scale



Close-up of the DS8150 display, showing the 1" tall LCD display and the easy to use controls.

FEATURES:

- EMR Ready
- BMI Calculator Function
- Armrests flip up to aid in patient entry and exit
- Hinged footrests for comfort and easier patient transfer
- Zero/Units/Reweigh/BMI
- USB and RS-232 Output
- Remaining battery life indicator
- Powered by six AA batteries, auto-shut-off timer maximizes battery life.
- Hold feature to lock patient's weight on the display.

MAXIMUM CAPACITY

550 lb x 0.2 lb (250 kg x 0.1 kg)

ACCURACY

+/- 0.2 lb (0.1 kg)

DISPLAYED UNITS

lb or kg

POWER

AC Adapter (Included)
or Six AA Batteries (Included)

CONSTRUCTION

Powder Coated Steel
and Heavy Duty ABS Plastic

SEAT DIMENSIONS

21" W x 18" H

WARRANTY

 Two Years



Doran Scales, Inc. 1-800-264-4106
1315 Paramount Parkway Fax: 630-879-0073
Batavia, IL 60510 www.doranmedical.com