



HCPCS
A7032
A7034
A7035



NasalFit Deluxe EZ CPAP Mask

Item #	Product Description	UOM
100NDES	NasalFit Deluxe EZ CPAP Mask, Small	1/ea
100NDEM	NasalFit Deluxe EZ CPAP Mask, Medium	1/ea
100NDEL	NasalFit Deluxe EZ CPAP Mask, Large	1/ea

Accessories		
100NDEH	NasalFit Deluxe EZ Replacement Headgear	1/ea
100NDES-NH	NasalFit Deluxe EZ CPAP Mask without Headgear, Small	1/ea
100NDEM-NH	NasalFit Deluxe EZ CPAP Mask without Headgear, Medium	1/ea
100NDEL-NH	NasalFit Deluxe EZ CPAP Mask without Headgear, Large	1/ea
100NDES-CUSHION	NasalFit Deluxe EZ Replacement Cushion, Small	1/ea
100NDEM-CUSHION	NasalFit Deluxe EZ Replacement Cushion, Medium	1/ea
100NDEL-CUSHION	NasalFit Deluxe EZ Replacement Cushion, Large	1/ea

- The NasalFit Deluxe EZ CPAP Mask has a soft silicon cushion, comfortable forehead plate with silicon pad, and a 360° rotational elbow, which are designed to increase comfort, compliance, and mask stability.
- Adjustable headgear, easy to assemble cushions, and a lightweight streamlined design help to ensure proper fit and minimal leak during sleep.
- By providing a mask option without headgear, a more cost effective option is available based upon the reimbursement schedule.
- Latex Free.