



10105-2

## Platform Walker/Crutch Attachment

- 10105-1 For use with Adult and Junior Walkers and Aluminum Crutches, 1/cs
- 10105-2 For use with Adult and Junior Walkers and Aluminum Crutches, 2/cs
- 10116 Universal Platform Walker/Crutch Attachment, 1/cs (For use with all Drive walkers and Invacare® Paddle Release Walkers)
- 10105C Replacement Clamps for 10105-2, 1/ea
- 10116C Replacement Clamps for 10116, 1/ea
- 10105CB Crutch Attachment Bracket for use with 10105-1, 10105-2, (Must be used with 10105-1, 10105-2), 1/ea

## Bariatric Platform Walker/Crutch Attachment

Available  
**1 Pack**

- 10105HD-2 Bariatric, 500 lb. Weight Capacity, Universal Platform Walker and Aluminum Crutch Attachment, 2/cs For Use on Item # 10220-2, 10220HD and 10223-2, 2/cs



### FEATURES:

- Designed for those who are unable to grip the walker.
- Cradles forearm with soft vinyl padding for comfortable weight bearing.
- Large top bracket allows the platform attachment to be placed over the hand grip of the walker.
- Adjustable strap securely holds arm in place.
- Limited Lifetime Warranty.
- Attachment for Underarm Crutch for item #10105 must be ordered separately. (Cannot be used with Bariatric model)



10105HD-2



10116

### CAN BE USED WITH CRUTCH



10105CB

### CAN ATTACH TO WALKER BRACE OR WALKER HAND GRIP



10105CB

Attachment for Underarm Crutch for item #10105 must be ordered separately. (Cannot be used with Bariatric model)

### SPECIFICATIONS

Item #	Overall Width	Overall Depth	Height (highest setting)	Height (lowest setting)	Product Weight	Carton Shipping Weight	Weight Limit
10105-1	6"	18"	50.5"	35.5"	3 lbs.	4.5 lbs.	300 lbs.
10105-2	6"	18"	50.5"	35.5"	3 lbs.	8 lbs.	300 lbs.
10116	5.5"	18"	50"	35.5"	4.5 lbs.	11 lbs.	300 lbs.
10105HD-2	6"	18"	52"	40"	4.5 lbs.	9 lbs.	500 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.