



Deluxe Folding Walker, One Button

10202-2 Adult, 2/cs
10203-2 Junior, 2/cs

FEATURES:

- Single-button release allows quick, one-hand folding.
- Cross brace is located near hand grips to allow for full stride.
- Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength while remaining lightweight.
- Limited Lifetime Warranty.



Deluxe Folding Walker, One Button with 5" Wheels

10212-2 Adult, 5" Wheels, 2/cs
10213-2* Junior, 5" Wheels, 2/cs

FEATURES:

- Same features as 10202-2/10203-2.
- Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength while remaining lightweight.
- Comes with wheels. For replacement wheels, refer to item # 10109 for 5" wheels.
- Includes rear glide caps, (item # 10107) and glide covers (item # 10107C) allowing use on all surfaces.
- Limited Lifetime Warranty.

SPECIFICATIONS

Item #	Width (inside hand grip)	Width (inside back legs)	Overall Width	Depth (opened at base)	Height (highest setting)	Height (lowest setting)	Product Weight	Carton Shipping Weight	Weight Limit
10202-2	17"	20.5"	23"	18.5"	39"	32"	6.0 lbs.	14 lbs.	300 lbs.
10203-2	17"	20.5"	23"	18.5"	32"	25"	5.8 lbs.	13.6 lbs.	300 lbs.
10212-2	17"	20.5"	25.5"	20.5"	39"	32"	7 lbs.	16 lbs.	300 lbs.
10213-2	17"	20.5"	25.5"	20.5"	32"	25"	6.8 lbs.	15.6 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

walkers

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM