

BARIATRIC
500 LB
WEIGHT CAPACITY



THIS ITEM IS AVAILABLE IN A 1-PACK!

EASY HEIGHT ADJUSTMENT WITH SELF THREADING KNOB!



10215RD-1
SEE PAGE 52 & 53 FOR REPLACEMENT SEATS.

Go-Lite Bariatric Steel Rollator, Padded Seat, 8" Casters with Loop Locks

10215BL-1 Blue, (Has replaced 685 series and 466), 1/cs
10215RD-1 Red, (Has replaced 685 series), 1/cs

FEATURES:


- Large 8" casters are ideal for indoor and outdoor use.
- Comes with large 21" (W) x 9" (D) x 6" (H) basket that can be mounted under seat (item # 10215BSK).
- Tool Free removable padded back rest for comfort.
- Strong Steel reinforced frame supports 500 lbs.
- Soft padded oversized seat (18" (W) x 13" (D)).
- Special loop lock made of internal aluminum casting operates easily and ensures safety.
- Handles are adjustable in height.
- Limited Lifetime Warranty.
Lifetime Warranty on brake cable.

Adjustable Width Winnie Walker, Bariatric EXP Steel Rollator, 8" Casters with Loop Locks

466* Light Blue & Chrome, 1/cs

FEATURES:

- Large 8" casters are ideal for indoor and outdoor use.
- 4" adjustable width allows entry into narrow doorways.
- Strong Steel reinforced frame supports 400 lbs.
- Special loop lock made of internal aluminum casting operates easily and ensures safety.
- Soft padded oversized seat.
- Handles are adjustable in height.
- Comes with oversized 14" x 12" x 7" wire basket (item # 828).
- Limited Lifetime Warranty.
Lifetime Warranty on brake cable.



PARTS DIAGRAMS ON PAGE 67

BARIATRIC
400 LB
WEIGHT CAPACITY

EASY TO REACH KNOB ADJUSTS ROLLATOR WIDTH 4" TO ALLOW ENTRY INTO NARROW DOORWAYS.



466*
SEE PAGE 52 & 53 FOR REPLACEMENT SEATS.

ADJUST FRAME WIDTH UP TO 4" WHEN USING LEVER



SEE ROLLATOR ACCESSORIES ON PAGES 50-57.
ROLLATOR PARTS DIAGRAMS ON PAGE 58-67.
ROLLATOR COLOR PAINT SWATCHES ON PAGE 68.

SPECIFICATIONS

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

Item #	Handle Height (highest)	Handle Height (lowest)	Unit Width	Unit Length	Seat (depth)	Seat (width)	Seat (height)	Product Weight	Carton Shipping Weight	Weight Limit
GO-LITE BARIATRIC	39"	34"	30"	29"	13"	18"	22"	26 lbs.	30.5 lbs.	500 lbs.
EXP/466	38"	34"	25" - 29"	30"	10"	20"	25"	28 lbs.	32.5 lbs.	400 lbs.

rollators

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM