



## "C" MODELS Glider Walker

COME WITH  
2 COLORS!



THIS ITEM  
IS AVAILABLE  
IN A 1-PACK!

10219C-2	Combo Pack, Standard, Includes 1 Flame Red and 1 Flame Blue, 2/cs
10219FRD-1	Standard, Flame Red, 1/cs
10219FBL-1	Standard, Flame Blue, 1/cs
10222C-2	Combo Pack, Junior, Includes 1 Flame Red and 1 Flame Blue, 2/cs
10222FRD-1	Junior, Flame Red, 1/cs
10222FBL-1	Junior, Flame Blue, 1/cs
10221C-2	Combo Pack, Pediatric, Includes 1 Flame Red and 1 Flame Blue, 2/cs
10221FRD-1	Pediatric, Flame Red, 1/cs
10221FBL-1	Pediatric, Flame Blue, 1/cs

### FEATURES:

- Aluminum frame with choice of colored cross bar.
- Comes with adjustable handles and legs that make this walker extremely versatile to accommodate virtually any physical requirements.
- Flat-free tires allow for both indoor and outdoor use.
- Wide base provides additional stability.
- Straight handle on standard and junior model; angled handle on pediatric model.
- Limited Lifetime Warranty.



## GLIDER WALKER ACCESSORIES

### Platform Attachment

10219PA	Standard, 1 pr/bx
10221PA	Pediatric, 1 pr/bx
10222PA	Junior, 1 pr/bx

### SPECIFICATIONS

Item #	Width (inside hand grip)	Width (inside back legs)	Depth (opened at base)	Height (highest setting)	Height (lowest setting)	Product Weight	Carton Shipping Weight	Weight Limit
10219FRD/FBL/C-1/2	Adjustable	28"	25"	37"	29"	8 lbs.	11/19 lbs.	175 lbs.
10221FRD/FBL/C-1/2	Adjustable	22"	20"	21"	12"	5.5 lbs.	8.5/15 lbs.	175 lbs.
10222FRD/FBL/C-1/2	Adjustable	24"	25"	33"	22"	7 lbs.	10/18 lbs.	175 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.

walkers  
To order call toll free: 877.224.0946

drive  
WWW.DRIVEMEDICAL.COM