



10224-4



10226-4

SHOWN WITH
5" WHEELS

Deluxe, Trigger Release Folding Walker

10224-4

Adult, Flame Blue Crossbrace,
(Has replaced 10224-2), 4/cs

10225-2*


Junior, Flame Blue Crossbrace, 2/cs

10225-4

Junior, Flame Blue Crossbrace, 4/cs

FEATURES:

- Trigger release feature allows the user to close the walker without removing their hands from the handle and is ideal for individuals with limited finger dexterity. Trigger release can be reached with hands of all sizes.
- Deluxe soft ribbed contoured hand grip for added comfort.
- Durable, Composite Trigger Release resists cracking and breaking.
- Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength while remaining lightweight.
- Limited Lifetime Warranty.
Patent # US6,886,575 B2



Flame Blue
Crossbrace



**NEW & IMPROVED
EASY "PULL UP
TRIGGER RELEASE"
FOR OPENING &
CLOSING**

Deluxe, Trigger Release Folding Walker with 5" Wheels

10226-4

Adult, Flame Blue with 5" Wheels, 4/cs

10227-2*


Junior, Flame Blue with 5" Wheels, 2/cs

10227-4

Junior, Flame Blue with 5" Wheels, 4/cs

FEATURES:

- Same features as 10224-4/10225-4.
- Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength while remaining lightweight.
- Replacement wheels 5" (10109).
- Includes rear glide caps, (item # 10107) and glide covers (item # 10107C) allowing use on all surfaces.
- Limited Lifetime Warranty.



Flame Blue
Crossbrace



**NEW & IMPROVED
EASY "PULL UP
TRIGGER RELEASE"
FOR OPENING &
CLOSING**

SPECIFICATIONS

Item #	Width (inside hand grip)	Width (inside back legs)	Overall Width	Depth (opened at base)	Height (highest setting)	Height (lowest setting)	Product Weight	Carton Shipping Weight	Weight Limit
10224-4	18"	21"	24"	18"	39"	32"	6.5 lbs.	28.5 lbs.	350 lbs.
10225-2/4	18"	21"	24"	17"	32.5"	25.5"	6.2 lbs.	14.8/27.3 lbs.	350 lbs.
10226-4	18"	21"	27"	21.5"	39.5"	32.5"	7.5 lbs.	28.5 lbs.	350 lbs.
10227-2/4	18"	20.5"	25.5"	19"	33.25"	26.5"	6.5 lbs.	15/28.5 lbs.	350 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.