

REPLACEMENT ACCESSORIES

SEAT AND BACKREST COVERS

460900252	Cover-Set Comfort, White, 1/bx
460900312	Cover-Set Comfort, Blue, 1/bx
460900432	Cover-Set Comfort, Grey, 1/bx

SPARE PARTS


460900403	Charger for Bellavita, 1/bx
415252	Plug + Washer for Cover, White, 1/ea
415312	Plug + Washer for Cover, Blue, 1/ea
415432	Plug + Washer for Cover, Grey, 1/ea
460900213	Suction Cup with clip, 2.95" (75 mm) Grey, 1/ea
460900220	Set of Suction Cups with Clips, 4 x 2.95" (75 mm) Grey, 4/pk
460900600	Hand Control Including Storage Battery, 1/ea

ACCESSORIES

101385	Rotating & Transfer Aid, White, 1/cs
102385	Rotating Aid, White, 1/cs
410200252	Headrest, Padded, White, 1/ea
410200312	Headrest, Padded, Blue, 1/ea
410200432	Headrest, Padded, Grey, 1/ea
460900205	Side Flap, 1/ea



bellavita



HIGHEST SEAT POSITION

LIGHTEST BATHLIFT
AVAILABLE

RECLINING BACKREST
LOWEST BATHLIFT AVAILABLE

SPECIFICATIONS:

Overall Length (with backrest in a reclined position)	32.6"
Overall Width (seat incl. side flaps)	27.5"
Width of Seat	13.7"
Depth of seat	16.5"
Overall Length of Seat	22.4"
Sitting Height Min.	2.3"
Sitting Height Max.	18.8"
Height of Backrest	25.1"
Width of Backrest	12.5" - 14.1"
Backrest Reclines	50°
Width of Baseplate	12.2"
Length of Baseplate	22"
Load Capacity Max.	300 lbs.
Total Weight (excl. hand control)	20.5 lbs.
Shipping Weight	25.3 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

bath safety

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM