



**THIS ITEM
IS AVAILABLE
IN A 1-PACK!**

Side Walker

10240-1 Adult, (Hemi-Walker), 1/cs

10240-2* Adult, (Hemi-Walker), 2/cs

FEATURES:

- Designed for individuals with the use of only one hand or arm.
- Lighter than a walker and more stable than a cane.
- Height adjustment 29.5" - 37".
- Folds easily with one hand.
- Ideal for users who do not require a walker but need a wider base for support.
- Limited Lifetime Warranty.



Bariatric Side Walker

10242 Adult, (Hemi-Walker), 1/cs

FEATURES:

- **Manufactured with reinforced tubing.**
- Designed for individuals with the use of only one hand or arm.
- Lighter than a walker and more stable than a cane.
- Height adjustment 30" - 36".
- Folds easily with one hand.
- Ideal for users who do not require a walker but need a wider base for support.
- Limited Lifetime Warranty.

SPECIFICATIONS

Item #	Overall Width	Overall Depth	Height (highest setting)	Height (lowest setting)	Product Weight	Carton Shipping Weight	Weight Limit
10240-1/2	17.5"	15.5"	37"	29.5"	3.5 lbs.	5.5/9 lbs.	300 lbs.
10242	16.5"	16"	36"	30"	10 lbs.	13 lbs.	500 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.