



New Closing Device Makes it Easy to Engage or Disengage.



LiftWalker™ with Retractable Stand Assist Bars and 5" Wheels

Folding Walker, Two Button with 5" Wheels

Item #	Product Description	UOM
10277LW	Two Button, Adult, with Retractable Assist Bars	1/cs

- LiftWalker™ will provide convenient, safe, assisted independence anywhere, anytime it's needed.
- For use by individuals who need assistance with sitting/standing on their own.
- Safer than being pulled by an attendant or caregiver.
- Reduces risk of injury to an attendant or caregiver.
- The retractable stand assist bars are easy to deploy. (Patented)
- The LiftWalker™ makes it easy to stand because the individual uses biceps, chest, back, stomach and leg muscles all at once. (Figure A - C)
- Easy to lock/unlock foot brake.

specifications	
Width	17" inside hand grip 21" inside back legs
Overall Width	26"
Depth (Opened at Base)	21"
Height	32.5" - 38.5"
Weight	11 lbs.
Carton Shipping Weight	13 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

Item #	Product Description	UOM
10244-1	Adult	1/cs
10244-4	Adult	4/cs

- Easy push-button release mechanisms.
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing.
- Sturdy 1" diameter construction ensures maximum strength.
- Newly designed rear glide cap allows for walker to slide easily and smoothly over most surfaces.
- Vinyl-contoured hand grip.
- Comes with wheels. (Replacement wheels item # 10109).
- Limited Lifetime Warranty.

specifications	
Width	17" inside hand grip 20" inside back legs
Overall Width	24"
Depth (Opened at Base) (10244-1 10244-4)	17.5"
Height (10244-1 10244-4)	32"-39"
Weight (10244-1 10244-4)	8 lbs.
Carton Shipping Weight	9 lbs./33 lbs.
Weight Limit	350 lbs.
Warranty	Limited Lifetime