



10406

**BARIATRIC**  
**500 LB**  
WEIGHT CAPACITY

## Bariatric Steel Crutch with Accessories

10406	Adult, 1pr/cs (500 lb. weight capacity)
10408	Tall Adult, 1pr/cs (500 lb. weight capacity)
10418*	Extra Tall Adult, 1pr/cs (550 lb. weight capacity)

### FEATURES:

- Durable, steel construction.
- Push-buttons easily adjust height in 1" increments.
- Comfortable, contoured underarm design.
- Large, non-skid jumbo size vinyl contoured tips.
- Limited Lifetime Warranty.



10268-1

COMPARE TO:

**drive** EZ Accessories®  
10268-1 EZ0015



**THIS ITEM  
IS AVAILABLE  
IN A 1-PACK!**

## Universal Cane/Crutch Pouch

10268*	3/cs
10268-1	1/cs

### FEATURES:

- Easily attaches to cane or crutch with adjustable Velcro-type closures.
- Stores personal items without interfering with mobility.
- Made of durable, easy to clean nylon.
- Size: 10" x 5" x 1".
- Warranty: 90 days.

### SPECIFICATIONS

Item #	Height (highest setting) (armrest to floor)	Height (lowest setting) (armrest to floor)	Patient Height (highest)	Patient Height (lowest)	Weight Limit
10406	55"	47"	5'10"	5'2"	500 lbs.
10408	63"	55"	6'6"	5'10"	500 lbs.
10418	66"	57"	6'8"	5'10"	550 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.

canes & crutches

To order call toll free: 877.224.0946

**drive**  
WWW.DRIVEMEDICAL.COM